

# Corning® Gosselin™ Square HDPE Storage Bottles

High density polyethylene bottles with wide neck

CORNING

- ▶ HDPE storage bottles are easier to handle, require less space (13% to 20%) on the shelf, and are ideal for sampling and storage
- ▶ HDPE resin is lightweight and break-resistant
- ▶ Bottles are leak-proof

## Key Features

- ▶ Good for storage of liquid, semi-liquid, and solid forms
- ▶ Available in six sizes: 150 mL, 250 mL, 500 mL, 1L, 2L, and 2.5L
- ▶ Large neck opening provides easier access for liquid transfer
- ▶ Graduated (except 150 mL)
- ▶ Sterile
- ▶ Raw materials are food contact compliant



## Ordering Information

### Corning® Gosselin™ Square HDPE Storage Bottles

Cat. No.	Description	Qty/Cs
LC150A-05	150 mL, 37 mm red screw cap, sterile	300
LC250-06	250 mL, graduated, 37 mm red screw cap, sterile	210
LC500-11	500 mL, graduated, 58 mm red screw cap, sterile	175
LC1000-13	1L, graduated, 58 mm red screw cap, sterile	90
LC2000-03	2L, graduated, 58 mm red screw cap, sterile	50
LC2500A-07	2.5L, graduated, 77 mm white screw cap, sterile	35

For more specific information on claims, visit the Certificates page at [www.corning.com/lifesciences](http://www.corning.com/lifesciences).

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

# CORNING

**Corning Incorporated**  
*Life Sciences*

836 North St.  
Building 300, Suite 3401  
Tewksbury, MA 01876  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

**ASIA/PACIFIC**  
**Australia/New Zealand**  
t 61 427286832

**China**  
t 86 21 3338 4338  
f 86 21 3338 4300

**India**  
t 91 124 4604000  
f 91 124 4604099

**Japan**  
t 81 3-3586 1996  
f 81 3-3586 1291

**Korea**  
t 82 2-796-9500  
f 82 2-796-9300

**Singapore**  
t 65 6572-9740  
f 65 6735-2913

**Taiwan**  
t 886 2-2716-0338  
f 886 2-2516-7500

**EUROPE**  
CSEurope@corning.com

**France**  
t 0800 916 882  
f 0800 918 636

**Germany**  
t 0800 101 1153  
f 0800 101 2427

**The Netherlands**  
t 020 655 79 28  
f 020 659 76 73

**United Kingdom**  
t 0800 376 8660  
f 0800 279 1117

**All Other European Countries**  
t +31 (0) 206 59 60 51  
f +31 (0) 206 59 76 73

**LATIN AMERICA**  
grupoLA@corning.com

**Brasil**  
t 55 (11) 3089-7400

**Mexico**  
t (52-81) 8158-8400