

Corning Incorporated  
Life Sciences

Registered  
ISO 9001:2008

## Product Description

**Catalog Number:** 1730-10

**Product Description:** Corning® Snap-Seal, Container, Sample, Opaque, Not Sterile, 10.0 OZ / 300 mL

**Component Materials:**  
Polypropylene (PP)

**Attributes:**

Capacity - 300 mL  $\pm$ 15 mL (5%)  
Product Dimensions - 63 mm OD x 112 mm Height

**Autoclaving:**

This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes. NOTE: Autoclaving disposable products IS NOT recommended.

**Heating / Cooling:**

This product can be heated to 140°C.

This product can be cooled to -196°C.

**Vacuum / Pressure:**

This product is not designed for and is not rated for vacuum or pressure applications.

**Sterilization:**

Shipped Not Sterile

**Lot Number Designation:**

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

**Revision Date:**

17 JUN 2013

Rev No: 1