

Corning® Lambda™ Plus Pipettors

Trusted design for your lab

CORNING



Single-channel pipettor



8-channel and 12-channel pipettors

Corning Lambda Plus Pipettors

Corning Lambda Plus pipettors have been engineered to provide the highest levels of accuracy, reproducibility, and comfort. Available in single-, 8- and 12-channel configurations, all Lambda Plus pipettors feature smooth plunger movement and very low pipetting forces. Volume can be set by turning either the thumb wheel or the plunger button, and a locking mechanism prevents accidental volume change.

The Lambda Plus pipettors have been designed to maintain their calibration, even after repeated use and autoclaving. A key is included with each pipettor for easy in-house calibration.

- ▶ Factory calibrated to strict ISO 8655 standards
- ▶ 3-year warranty
- ▶ Fully autoclavable and UV resistant
- ▶ Designed for universal fit with all common brands of pipet tips
- ▶ Adjustable tip ejector

Corning Lambda Plus Starter Kit

The Lambda Plus starter kit contains four Lambda Plus single-channel pipettors conveniently packaged with useful accessories (manual, calibration key, color-coded rings, tube of piston grease, and identification stickers). The kit includes:

- ▶ 4 Corning Lambda Plus single-channel pipettors:
 - 0.5 to 10 μL (Cat. No. 4071)
 - 2 to 20 μL (Cat. No. 4072)
 - 20 to 200 μL (Cat. No. 4074)
 - 100 to 1000 μL (Cat. No. 4075)
- ▶ Universal carousel rack for six single-channel pipettors (Cat. No. 4078)
- ▶ Corning DeckWorks pipet tips: Cat. Nos. 4135 (10 μL), 4138 (200 μL), and 4140 (1000 μL)



Starter kit

Ordering Information

Corning® Lambda™ Plus Single-channel, 8-channel, and 12-channel Pipettors

Cat. No.	Volume Range (µL)	Channels	Accuracy (%)	Precision (%)	Non-filtered Tips (µL)	Qty/Cs
4070	0.1 - 2.0	1	±40.0 to ±1.5	≤12.0 to ≤0.7	10	1
4071	0.5 - 10	1	±4.0 to ±0.5	≤2.8 to ≤0.4	10	1
4072	2 - 20	1	±3.0 to ±0.8	≤1.5 to ≤0.3	200	1
4073	10 - 100	1	±1.6 to ±0.8	≤0.8 to ≤0.2	200	1
4074	20 - 200	1	±1.2 to ±0.6	≤0.6 to ≤0.2	200	1
4075	100 - 1,000	1	±1.6 to ±0.6	≤0.40 to ≤0.15	1,000	1
4080	0.5 - 10	8	±10.0 to ±2.0	≤8.0 to ≤1.2	10	1
4081	5 - 50	8	±4.0 to ±1.6	≤2.5 to ≤0.6	200	1
4082	20 - 200	8	±3.0 to ±1.0	≤1.5 to ≤0.6	200/300	1
4083	50 - 300	8	±1.6 to ±1.0	≤1.5 to ≤0.6	300	1
4084	0.5 - 10	12	±10.0 to ±2.0	≤8.0 to ≤1.2	10	1
4085	5 - 50	12	±4.0 to ±1.6	≤2.5 to ≤0.6	200	1
4086	20 - 200	12	±3.0 to ±1.0	≤1.5 to ≤0.6	200/300	1
4087	50 - 300	12	±1.6 to ±1.0	≤1.5 to ≤0.6	300	1



Corning Lambda Plus Starter Kit

Cat. No.	Description	Qty/Cs
4069	Corning Lambda Plus starter kit Includes: 4 single-channel pipettors: 0.5 to 10 µL (Cat. No. 4071), 2 to 20 µL (Cat. No. 4072), 20 to 200 µL (Cat. No. 4074), 100 to 1000 µL (Cat. No. 4075); Universal carousel rack for six single-channel pipettors (Cat. No. 4078); Corning DeckWorks pipet tips: 10 µL (Cat. No. 4135), 200 µL (Cat. No. 4138), 1000 µL (Cat. No. 4140)	1 kit

Accessories

4078	Universal carousel rack for six Corning Lambda Plus single-channel pipettors	1
4079	Universal linear rack for six Corning Lambda Plus single- and multi-channel pipettors	1
4088	Pipettor stand for 8- and 12-channel Corning pipettors	1
6066	Universal linear stand for four pipettors, orange (four single-channel, maximum of two 8-channel or 12-channel pipettors)	1
6067	Universal linear stand for four pipettors, transparent (four single-channel, maximum of two 8-channel or 12-channel pipettors)	1
6068	Universal linear stand for six pipettors, transparent (six single-channel, maximum of three 8-channel or 12-channel pipettors)	1

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

CORNING

Corning Incorporated
Life Sciences

www.corning.com/lifesciences

NORTH AMERICA

t 800.492.1110
t 978.442.2200

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832

Chinese Mainland
t 86 21 3338 4338

India

t 91 124 4604000

Japan

t 81 3-3586 1996

Korea

t 82 2-796-9500

Singapore

t 65 6572-9740

Taiwan

t 886 2-2716-0338

EUROPE

CSEurope@corning.com

France

t 0800 916 882

Germany

t 0800 101 1153

The Netherlands

t 020 655 79 28

United Kingdom

t 0800 376 8660

All Other European Countries

t +31 (0) 206 59 60 51

LATIN AMERICA

grupoLA@corning.com

Brazil

t 55 (11) 3089-7400

Mexico

t (52-81) 8158-8400