

Corning Incorporated  
Life Sciences

Registered  
ISO 9001:2008

## Product Description

**Catalog Number:** 2949-75X25

**Product Description:** Corning® Microscope Slide, Frosted Two Sides, One End, 75 x 25 mm

**Component Materials:**

Type II - Soda-Lime Glass, per ASTM E-438 (Standard Specification for Glasses in Laboratory Apparatus)

**Attributes:**

Product Dimensions - 75 mm x 25 mm x 0.90 - 1.10 mm Thick

**Autoclaving:**

This product is made of Soda Lime glass and IS NOT designed to be autoclaved.

**Heating / Cooling:**

This product is designed for use in general laboratory microscopy. It is not intended for heating or cooling.

**Vacuum / Pressure:**

This product is not designed for and is not rated for vacuum or pressure applications.

**Sterilization:**

Shipped Not Sterile

**Lot Number Designation:**

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 1