



The following features are present in the vector based on the nucleotide sequence.

Thymidine kinase (TK) promoter .....	26-777
Green <i>Renilla</i> luciferase gene .....	791-1725
BGH poly(A) signal.....	1734-1859
SV40 origin/promoter.....	1860-2424
Puromycin resistant gene .....	2425-3024
SV40 poly(A) signal .....	3186-3219
beta-lactamase (Amp <sup>R</sup> ) gene .....	3328-4188
pUC replication origin (pUC Ori).....	4367-5171
Transcriptional terminator (Ter) .....	5172-5779
Lac operator 1 (Lac O1) .....	5780-5800
Transcriptional pause site (TPS).....	5933-6004