

**Corning Incorporated  
Life Sciences**

**Registered  
ISO 9001:2008**

## Product Description

**Catalog Number:** 1395-32HTC

**Product Description:** Corning® Screw Cap, High Temperature, PBT with PTFE Liner, GL-32

**Component Materials:**

Cap - Polybutylene terephthalate (PBT), Color: Red  
Liner - Polytetrafluoroethylene (PTFE) Faced Silicone

**Attributes:**

Product Dimensions - 42 mm OD x 26 mm High

**Autoclaving:**

This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes.

**Heating / Cooling:**

The normal use temperature range for this product is from -50°C to +125°C. The upper temperature limit for the plastic components of this product is 180°C. This product is not designed for direct contact heating.

**Vacuum / Pressure:**

This product is a component part and not specifically rated for vacuum or pressure.

**Sterilization:**

Shipped Not Sterile

**TSE/BSE:**

Materials of animal origin are not used in the manufacture of this product.

**Lot Number Designation:**

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

**Revision Date:**

14 JULY 2014

Rev No: 4