

Split-Type Flake Ice Maker



FIM300



FIM500

Applicable area

Various types of experiments, hospitals, schools, laboratories, medical, chemical applications, supermarkets, hotels and other seafood products stores etc.

Technical Parameters

Model	FIM300	FIM500
Ice Making Capacity(kg/24h)	300	500
Ice Storage Capacity(kg)	150	200
Ice Shape	Small particle of irregular snow ice	
Input Current(w)	1360	2100
Way of Condensation	Air cooling	
Ambient Temperature	10~35℃	
Water Temperature	2~15℃	
Tank Shell	Stainless steel	
Refrigerant	SECCOP /R404a	Embraco/R404a
Standard Accessories	Water inlet pipe*1, Drainage pipe*1, Ice spoon*1	
Optional Accessories	Filter	
Power Supply	AC 220V±10%,50/60HZ ; 110V±10%, 60HZ	
Net Weight(kg)	95	125
Gross Weight(kg)	125	155
External Size (W*D*H)mm	Head Size	600*550*570
	Ice Tank Size	570*830*1020
Package Size (W*D*H)mm	Head Size	640*860*890
	Ice Tank Size	640*860*990
External Size (W*D*H)mm	Head Size	800*760*740
	Ice Tank Size	770*765*11130
Package Size (W*D*H)mm	Head Size	700*850*995
	Ice Tank Size	880*870*1050