

Catalog Number: 431479

Product Description: Corning® Spin-X® UF Concentrator, 500 uL, 30k MWCO,
Non-Sterile

Component Materials:

- Inner Tube - Virgin Polycarbonate, meets *USP, Class VI* requirements for plastic containers and closures.
- Outer Tube - Virgin Polypropylene.
- Membrane - Polyethersulfone (PES).

Product Dimensions:

Tube Assembly Height - 1.97 inches Tube Diameter at Base - 0.43 inches

Maximum RCF:

Swing out (bucket) Not recommended
Fixed Angle (40-45°) 12,000xg

Membrane Retention:

The product has been tested for retention using a standard molecular weight marker solution.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Integrity Test: Pressure hold test method.
Visual Attributes: Visual examination of product.
Packaging: Visual examination of accurate labeling and correct product configuration.

Revision Date:

03-14-19

Rev No: 2