

**Corning Incorporated  
Life Sciences**

**Registered  
ISO 9001:2008**

## Product Description

**Catalog Number:** 431121

**Product Description:** Corning® cryogenic storage box for 1-2 mL cryogenic vials, holds 100 internally threaded vials

**Component Materials:**

Box - Virgin Polycarbonate  
Lid - Virgin Polycarbonate

**Capacity:**

Maximum capacity is 100 internally threaded vials.

**Temperature Range:**

Can withstand temperatures of - 196°C to +121°C. It is strongly recommended that you test under actual conditions before use. The unit is autoclavable.

**Bovine Spongiform Encephalopathy and Transmissible Spongiform Encephalopathy:**

This product complies with the latest revision of EMEA/410/01 "Note for Guidance on minimizing the risk of transmitting animal spongiform encephalopathy agents via human veterinary medicine products" by virtue of all bovine derived material having been processed per specific condition of section 6.4 of EMEA/410/01.

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes: Visual examination of the product.

Packaging: Inspection for accurate labeling and correct product configuration.

**Lot Number Designation:**

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Run number for that year.

Rev No: 3