



## Mini Centrifuge

### Features:

- \* Timing function.
- \* Rotor with 8 positions.
- \* Completely transparent cover.
- \* Humanization design, switch control for cover.
- \* Auto stop working when cover opens during running and restart running after cover closed.
- \* Convex type switch, easy to use.
- \* Low noise, stable running.
- \* Long life, reliable motor.
- \* Unique rotor clip design, easy to install and replace.
- \* Large LED display, indicate speed and set time.



### Technical Parameters:

Model	Mini-4K	Mini-6K	Mini-7K	Mini-10K	Mini-10K+
Max. Speed	4000rpm	6000rpm	7000rpm	10000rpm	3000-10000rpm
Max. RCF	860xg	2000xg	2600xg	5400xg	500-5400xg
Rotor	8*2.0ml/1.5ml/0.5ml/0.2ml 2*8*0.2ml PCR centrifugal tube		8*2.0ml/1.5ml/0.5ml/0.2ml		8*2.0ml/1.5ml/0.5ml/0.2ml 2*8*0.2ml PCR centrifugal tube(speeds≤6000rpm)
Noise	≤45dB		≤50dB		
Power Supply	AC110/220V±10%, 50/60Hz				
External Size(W*D*H)mm	176*160*121				
Package Size(W*D*H)mm	260*200*150				
Net/Gross Weight(kg)	1.4/3				

## Mini High Speed Centrifuge

### Technical Parameters:

Model	Mini-15K
Capacity	12*1.5ml/2.0ml
Speed Range	800~15000rpm
Max. RCF	15100xg
Timing Range	15s~99min
Consumption	85W
Power Supply	AC110/220V±10%, 50/60Hz
External Size(W*D*H) mm	230*220*130
Package Size(W*D*H) mm	330*340*270
Net Weight/Gross Weight(kg)	3.8/6



## Three Speed Mini Centrifuge

### Technical Parameters:

Model	Mini-4610KII
Speed Range	4000rpm/6000rpm/10000rpm
Max. RCF	770xg/1700xg/4800xg
No.1 Angle Rotor	6*2.0ml/1.5ml/0.5ml/0.2ml
No.2 Angle Rotor	2*8*0.2ml
Power Supply	DC12V
External Size(W*D*H) mm	170*150*130
Package Size(W*D*H) mm	440*400*360
Gross Weight(kg)	6



## Single Speed Mini Centrifuge

### Technical Parameters:

Model	Mini-4KH	Mini-6KH	Mini-10KH
Max. Speed	4000rpm	6000rpm	10000rpm
Max. RCF	770xg	1700xg	4800xg
No.1 Angle Rotor	6*2.0ml/1.5ml/0.5ml/0.2ml		
No.2 Angle Rotor	2*8*0.2ml		
Power Supply	AC110/220V±10%, 50/60Hz		
External Size(W*D*H) mm	170*150*130		
Package Size(W*D*H) mm	440*400*360		
Gross Weight(kg)	6		

