

Product Description

Catalog Number: 3936

Product Description: Corning® 1536-well, 10µl, non-treated, solid black plate, round bottom,

Component Materials:

Plate - Virgin Polystyrene, meets *USP, Class VI* requirements for plastic containers and closures. Black concentrate

Product Dimensions:

Length of Plate	-	5.030 in.	Width of Square Well @ Top	-	.0665 in.
Width of Plate	-	3.365 in.	Diameter of Well @ bottom	-	.029 in.
Depth of Well	-	.197 in.	Height	-	.410 in.
Tolerances of Dimensions	-	+/- .010 in.	Recommended working volume per well is	-	up to 8 µl

Surface Characterization:

Surface is non-treated and hydrophobic.

Optical Characteristics:

The product is made of opaque black polystyrene walls to minimize well to well crosstalk and background fluorescence.

Binding Characteristics:

The product has been tested for the attribute of uniform well to well protein binding utilizing a fluorescent assay.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes: Visual examination of the product.

Packaging: Inspection for accurate labeling, and correct product configuration.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian Date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.