


















Format: 96 tip
 Tip Volume Range Selection: 250 µL - 1100 µL
 Material: Clear Polypropylene
 Certified RNase-/DNase-free, nonpyrogenic
 Available Options:
 • With or without aerosol-resistant filters
 • Nonsterile or sterilized



Tip Box Format: Non-ANSI hanging tip rack
 Material: Number 7 Polycarbonate

Axygen® Robotic Tips

Type	Max Volume	Compatible Head	Tip Box	Packaging	Standard Cat. No.
Non-Filtered	300 µL	BioRobot		Nonsterile Sterile Quantity: 96 tips/rack, 24 racks/case	QT-300-R QT-300-R-S
					
Filtered	250 µL	BioRobot		Sterile Quantity: 96 tips/rack, 24 racks/case	QTF-300-R-S
					
Conductive, Non-Filtered	300 µL	BioRobot		Nonsterile Sterile Quantity: 96 tips/rack, 24 racks/case	QT-300-CBK-R QT-300-CBK-R-S
					
Conductive, Filtered	250 µL	BioRobot		Sterile Quantity: 96 tips/rack, 24 racks/case	QTF-300-CBK-R-S
					
Non-Filtered	1100 µL	BioRobot		Nonsterile Sterile Quantity: 96 tips/rack, 16 racks/case	QT-1100-R QT-1100-R-S
					
Filtered	1000 µL	BioRobot		Sterile Quantity: 96 tips/rack, 16 racks/case	QTF-1100-R-S
					
Conductive, Non-Filtered	1100 µL	BioRobot		Nonsterile Sterile Quantity: 96 tips/rack, 16 racks/case	QT-1100-CBK-R QT-1100-CBK-R-S
					
Conductive, Filtered	1000 µL	BioRobot		Sterile Quantity: 96 tips/rack, 16 racks/case	QTF-1100-CBK-R-S
					

To order Corning products, call 800.492.1110 or visit www.corning.com/lifesciences for a list of distributors.