
eBioscience™ IHC Antigen Retrieval Solution – High pH (10X)

Catalog Number: 00-4956

For Research Use Only. Not for use in diagnostic procedures.

Product Information

Contents: eBioscience™ IHC Antigen Retrieval Solution – High pH (10X)



Catalog Number: 00-4956

Formulation: 10X: 100 mM Tris, 10 mM EDTA solution, pH 9.0



Temperature Limitation: Store at room temperature.



Batch Code: Refer to vial



Use By: Refer to vial

Description

The IHC Antigen Retrieval Solution- High pH (10X) is designed for use during the heat induced epitope retrieval (HIER) step prior to immunohistochemistry on formalin-fixed paraffin embedded tissue sections. The use of this buffer in combination with heat (often by microwave, water bath, or pressure cooker) has been shown to restore the antigenicity of proteins modified during the formalin fixation of tissue. This buffer is supplied as a 10X stock solution. Dilution of this buffer to 1X using reagent grade water results in a 10 mM Tris, 1 mM EDTA solution, pH 9.0.

Applications Reported

IHC Antigen Retrieval Solution – High pH (10X) has been reported for use in immunohistochemical staining of formalin-fixed paraffin embedded tissue sections.

Applications Tested

IHC Antigen Retrieval Solution – High pH (10X) has been tested immunohistochemistry on formalin-fixed paraffin embedded tissue.

References

Lyck L, Dalmau I, Chemnitz J, Finsen B, Schrøder HD. Immunohistochemical markers for quantitative studies of neurons and glia in human neocortex. J Histochem Cytochem. 2008 Mar;56(3):201-21

Shi SR, Cote RJ, Taylor CR. Antigen retrieval techniques: current perspectives. J Histochem Cytochem. 2001 Aug;49(8):931-7.

Shi SR, Key ME, Kalra KL. Antigen retrieval in formalin-fixed, paraffin-embedded tissues: an enhancement method for immunohistochemical staining based on microwave oven heating of tissue sections. J Histochem Cytochem. 1991 Jun;39(6):741-8

Related Products

00-4955 eBioscience™ IHC Antigen Retrieval Solution – Low pH (10X)

Not for further distribution without written consent.

Copyright © 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

Tel: 888.999.1371 or 858.642.2058 • Fax: 858.642.2046 • thermofisher.com/ebioscience •

info@ebioscience.com