

## **TECHNICAL SHEET**

PRODUCT NAME:	Anti-Giardia mAb cyst protein (clone GR7)
Code:	MT-16GR7
Clone:	GR7
lsotype:	lgG1
Physical State:	Liquid
Host Animal (source):	Mouse
Specificity:	GR7 is specific for Giardia detection.
Cross Reactivity:	No cross reactivity with Cryptosporidium parvum, Entamoeba histolytica, Entamoeba dispar.
Purification:	Anti-Giardia mAb cyst protein (clone GR7) has been purified in a multi-step procedure. ≥95% pure by SDS-PAGE.
Storage Conditions:	Storage at 2-8ºC.
Health & Safety Inform	This product contains sodium azide as preservative. Although the amount is 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.
Presentation:	PBS pH: 7.4, 0.09% Sodium Azide (NaN₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 11/01/2019





Pol. Industrial Río Gállego II, Calle J, № 1, 50840, San Mateo de Gállego, Zaragoza (SPAIN) Tel. +34 976520354 <u>certest@certest.es</u> <u>www.certest.es</u>