

Introducing The New Corning® Hinged Rack Pipet Tips



CORNING
Discovering Beyond Imagination

Life
Sciences

New!

Corning Universal Fit Hinged Rack Pipet Tips are designed to assure a reliable fit with all major brand pipettors. They are available in two sizes: 200 μL and 1,000 μL .

The hinged rack allows easy rack stacking and steam penetration during autoclaving. Racks can be reused. Simply remove the empty tip insert in the rack and replace it with a new Corning tip insert. The inserts in the chart on the back of this sheet are compatible with hinged racks.

- ▶ 200 μL racks have 96 tips per rack;
1,000 μL racks have 100 tips per rack
- ▶ Available sterile and nonsterile
- ▶ Tips are not graduated
- ▶ Certified RNase-/DNase-free and nonpyrogenic
- ▶ Replacement inserts for hinged racks are a great way to reduce waste and reuse racks.



Corning® Hinged Rack Pipet Tips

Corning Hinged Rack Pipet Tips Ordering Information

Cat. No.	Volume Range	Format	Sterile	Color	Racks and/or Inserts/ Carton	Tips/ Carton
4711	1-200 µL	96 Tip Hinged Rack	Yes	Yellow	10	960
4712	1-200 µL	96 Tip Hinged Rack	No	Yellow	10	960
4710	1-200 µL	96 Tip Insert for Hinged Rack	No	Yellow	10 Inserts	960
4714	100-1000 µL	100 Tip Hinged Rack	Yes	Blue	10	1,000
4713	100-1000 µL	100 Tip Hinged Rack	No	Blue	10	1,000
4715	100-1000 µL	100 Tip Insert for Hinged Rack	No	Blue	10 Inserts	1,000

Empty hinged racks are not sold separately.

For additional product or technical information, please visit www.corning.com/lifesciences or call 800.492.1110. Customers outside the United States, call +1.978.635.2200 or contact your local Corning sales office listed below.

Corning Incorporated Life Sciences

45 Nagog Park
Acton, MA 01720
t 800.492.1110
t 978.635.2200
f 978.635.2476

www.corning.com/lifesciences

Worldwide Support Offices

ASIA / PACIFIC

Australia
t 61 2-9416-0492
f 61 2-9416-0493

China
t 86 21-3222-4666
f 86 21-6288-1575

Hong Kong
t 852-2807-2723
f 852-2807-2152

India

t 91-124-235 7850
f 91-124-401 0207

Japan

t 81 (0) 3-3586 1996/1997
f 81 (0) 3-3586 1291/1292

Korea

t 82 2-796-9500
f 82 2-796-9300

Singapore

t 65 6733-6511
f 65 6861-2913

Taiwan

t 886 2-2716-0338
f 886 2-2716-0339

EUROPE

France

t 0800 916 882
f 0800 918 636

Germany

t 0800 101 1153
f 0800 101 2427

The Netherlands & All Other

European Countries
t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

United Kingdom

t 0800 376 8660
f 0800 279 1117

LATIN AMERICA

Brasil

t (55-11) 3089-7419
f (55-11) 3167-0700

Mexico

t (52-81) 8158-8400
f (52-81) 8313-8589