

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 6120P-30

Product Description: Corning® Funnel, Reusable Plastic, PP, Not Sterile, 30 mm

Component Materials:
Polypropylene (PP)

Attributes:
Product Dimensions - 30 mm Top ID, 7 - 5 mm Stem ID, 45 mm High

Autoclaving:
Autoclaving this product using a standard cycle of 121°C and 15 psi for 15 minutes may change the physical shape of the product and its original volumetric accuracy.

Heating / Cooling:
The upper temperature limit for this plastic material is 125°C.

This product IS NOT suitable for Dry Heat Sterilization as it will melt. This product MUST NOT be heated by any direct, surface contacting heat source.

A lower limit usage temperature has not been established for this product.

Vacuum / Pressure:
This product is not designed for and is not rated for vacuum or pressure applications.

Sterilization:
Shipped Not Sterile

Lot Number Designation:
Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Revision Date:
01 OCT 2012

Rev No: 1