

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 1395-25HTSC

Product Description: Corning® Screw Cap, High Temperature, PBT, Open Center, GL-25

Component Materials:

Cap - Polybutylene terephthalate (PBT), Color: Red

Attributes:

Product Dimensions - 33 mm OD x 23 mm High, 15 mm Center Hole Diameter

Autoclaving:

This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes.

Heating / Cooling:

The normal use temperature range for this product is from -50°C to +125°C. The upper temperature limit for the plastic components of this product is 180°C. This product is not designed for direct contact heating.

Vacuum / Pressure:

This product is a component part and not specifically rated for vacuum or pressure.

Sterilization:

Shipped Not Sterile

TSE/BSE:

Materials of animal origin are not used in the manufacture of this product.

Lot Number Designation:

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Revision Date:

14 JULY 2014

Rev No: 1