



00-2018 1 Liter  
Stabilizer Solution  
for HRP  
Lot No.

### STABILIZER SOLUTION FOR HRP CONJUGATE

**FORM:**

Liquid. Clear solution. The product is a Tris buffer, pH 7.5 containing protein to stabilize the enzyme conjugate and 0.1% Proclin (isothiazolin derivatives) as preservative.

**SUGGESTED USAGE:**

This product can be used as a conjugate diluent in peroxidase-based immunoassays. Dilute the enzyme conjugate at a minimum of 1:50 dilution to ensure the stabilizing benefit of this solution.

**STORAGE:**

Store at 2-8°C.

**For Research Use Only.**

**CAUTION: Not for human or animal therapeutic or diagnostic use.**

[www.invitrogen.com](http://www.invitrogen.com)

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288

E-mail: [techsupport@invitrogen.com](mailto:techsupport@invitrogen.com)

PI002018

(Rev 11/09) DCC-09-1726

**Important Licensing Information** - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, [www.invitrogen.com](http://www.invitrogen.com)). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.



00-2018 1 Liter  
Stabilizer Solution  
for HRP  
Lot No.

### STABILIZER SOLUTION FOR HRP CONJUGATE

**FORM:**

Liquid. Clear solution. The product is a Tris buffer, pH 7.5 containing protein to stabilize the enzyme conjugate and 0.1% Proclin (isothiazolin derivatives) as preservative.

**SUGGESTED USAGE:**

This product can be used as a conjugate diluent in peroxidase-based immunoassays. Dilute the enzyme conjugate at a minimum of 1:50 dilution to ensure the stabilizing benefit of this solution.

**STORAGE:**

Store at 2-8°C.

**For Research Use Only.**

**CAUTION: Not for human or animal therapeutic or diagnostic use.**

[www.invitrogen.com](http://www.invitrogen.com)

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288

E-mail: [techsupport@invitrogen.com](mailto:techsupport@invitrogen.com)

PI002018

(Rev 11/09) DCC-09-1726

**Important Licensing Information** - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, [www.invitrogen.com](http://www.invitrogen.com)). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.