



The following features are present in the vector based on the nucleotide sequence.

Multiple cloning site (MCS)	1-42
Green <i>Renilla</i> luciferase gene	49-984
BGH poly(A) signal	993-1118
SV40 origin/promoter	1119-1683
Puromycin resistant gene	1684-2283
SV40 poly(A) signal	2445-2478
beta-lactamase (Amp <sup>R</sup> ) gene	2587-3447
pUC replication origin (pUC Ori)	3626-4430
Transcriptional terminator (Ter)	4431-5038
Lac operator 1 (Lac O1)	5039-5059
Transcriptional pause site (TPS)	5192-5263