

Cylinder

Graduated cylinder: tall form

Capacity [ml]	Subdivision [ml]	Height [mm]	Outer diameter [mm]	Error limit [± ml]		
				BLAUBRAND®	BLAUBRAND® USP	SILBERBRAND
5	0.1	115	12	0.05	0.05	0.08
10	0.2	140	14	0.10	0.10	0.15
25	0.5	170	21	0.25	0.17	0.4
50	1	200	25	0.5	0.25	0.8
100	1	260	29	0.5	0.5	0.8
250	2	335	39	1.0	1.0	1.5
500	5	365	53	2.5	2.0	4
1000	10	465	65	5	3.0	8
2000	20	505	85	10	6.0	15

Graduated cylinder: low form

Capacity [ml]	Subdivision [ml]	Height [mm]	Outer diameter [mm]	Error limit [± ml]
SILBERBRAND				
10	1	90	21	0.3
25	1	115	25	0.5
50	2	145	29	1
100	2	165	39	1
250	5	195	54	2
500	10	250	65	5
1000	20	285	85	10
2000	50	340	105	20

Mixing cylinder: tall form

Capacity [ml]	Subdivision [ml]	Height [mm]	Outer diameter [mm]	Error limit [± ml]	
				BLAUBRAND®	SILBERBRAND
10	0.2	160	14	0.10	0.15
25	0.5	190	21	0.25	0.4
50	1	220	25	0.5	0.8
100	1	285	29	0.5	0.8
250	2	350	39	1.0	1.5
500	5	395	53	2.5	4
1000	10	500	65	5	8

