

## Touch Screen Constant-Temperature Incubator

### Features:

- \* Large LCD touch screen.
- \* D type: with inner glass door.  
G type: with viewing window.
- \* Microprocessor controller.
- \* Stainless steel chamber.



Double-layer Door (D type)



Glass Viewing Window Door (G type)

### Safety Functions:

- \* Over temperature alarm.
- \* Independent limit temperature alarm system.
- \* Leakage current and over-voltage protector.

### Standard Accessories:



Test port



Base stand



Universal casters



UV lamp

### Optional Accessories:

### Technical Parameters:

Model	BJPX-H54BK(D/G)	BJPX-H88BK(D/G)	BJPX-H160BK(D/G)
Capacity	54L	88L	160L
Temperature Control Method	PT100		
Temperature Range	RT+5~85°C		
Temperature Stability	±0.5°C (at 37°C)		
Temperature Uniformity	±1°C (at 37°C)		
Ambient Temperature	5~30°C, recommend 25°C±2°C		
Timing Range	0~999h		
Shelves	2pcs/adjustable		
Power Consumption	345W		495W
Power Supply	AC220V±10%, 50/60Hz(Standard); AC110±10%, 60Hz(Optional)		
Internal Size(W*D*H)mm	395*346*395	500*364*500	558*518*600
External Size(W*D*H)mm	D type: 536*588*681 G type: 522*626*587	D type: 650*573*793 G type: 627*621*692	D type: 718*732*890 G type: 685*779*790
Package Size(W*D*H)mm	D type: 606*648*818 G type: 592*676*724	D type: 720*630*930 G type: 697*671*829	D type: 788*792*1027 G type: 755*829*927
Gross Weight(kg)	D type: 51 G type: 44	D type: 65 G type: 51	D type: 89 G type: 69