

Loading-side Panel:



7.0" SIEMENS color display touch screen  
CPU226 SIEMENS PLC as central processor

- Display Information:**
- \* Jacket pressure P1
  - \* Chamber pressure P2
  - \* Chamber temperature T2
  - \* Door operation
  - \* Printer ON/OFF
  - \* Operator ID and password
  - \* Alarm display and record
  - \* Program selection
  - \* System parameter setting
  - \* Manual operation
  - \* Cycle process and values in real time
- A variety of program:**
- \* Fabric 121°C and 134°C
  - \* Instrument 121°C and 134°C
  - \* Liquid 121°C
  - \* Rubber 121°C
  - \* Flash 134°C
  - \* Gravity 121°C
  - \* Automatic leak test
  - \* User-defined
  - \* Demo

Unloading-side Panel:



- Humanistic indicating system:**
- \* Back door closing indicating lamp
  - \* Back door closed lamp
  - \* Cycle running lamp
  - \* Cycle end indicating lamp
  - \* Breakdown indicating lamp
  - \* Power lamp

Sterilization consumables:



Model	Meaning
350	Chamber volume(L)
S/D	Door: S-double door, D-single door
No marks/D	Steam supply: no marks-outer steam supply(optional); D-built-in steam generator
A/B	Installation: A-Shelf loading; B-Loading by cart and trolley

Technical Parameters:

Model	Volume	Chamber Size (W*H*L)	External Size (W*H*L)	Net Weight	Power Consumption	Power Supply			
BKQ-140D-A	146L	420*420*830mm	1290*1680*1050mm	425 Kg	2.5KW	380V, 50Hz			
BKQ-140DD-A			525 Kg	18KW					
BKQ-140S-A			510 Kg	2.5KW					
BKQ-140SD-A			610 Kg	18KW					
BKQ-240D-A	241L	600*600*740mm	1370*1880*950mm	580 Kg	2.5KW	380V, 50Hz			
BKQ-240DD-A			660 Kg	24KW					
BKQ-240S-A			660 Kg	2.5KW					
BKQ-240SD-A			720 Kg	24KW					
BKQ-350D-A	352L	600*600*980mm	1215*1880*1190mm	780 Kg	2.5KW	380V, 50Hz			
BKQ-350DD-A			880 Kg	24KW					
BKQ-350S-A			860 Kg	2.5KW	380V, 50Hz				
BKQ-350SD-A			960 Kg	24KW					
BKQ-350D-B			880 Kg	2.5KW					
BKQ-350DD-B			1000 Kg	24KW					
BKQ-350S-B			960 Kg	2.5KW					
BKQ-350SD-B			1080 Kg	24KW					
BKQ-650D-B			655L	610*910*1180mm		1310*2070*1370mm	1170 Kg	4KW	380V, 50Hz
BKQ-650DD-B						1270 Kg	46KW		
BKQ-650S-B	1270 Kg	4KW							
BKQ-650SD-B	1370 Kg	46KW							
BKQ-810D-B	810L	610*910*1460mm	1310*2070*1650mm	1260 Kg	4KW	380V, 50Hz			
BKQ-810DD-B			1360 Kg	61KW					
BKQ-810S-B			1360 Kg	4KW					
BKQ-810SD-B			1460 Kg	61Kw					
BKQ-990D-B	999L	610*910*1800mm	1310*2070*1990mm	1300 Kg	4KW	380V, 50Hz			
BKQ-990DD-B			1400 Kg	61KW					
BKQ-990S-B			1400 Kg	4KW					
BKQ-990SD-B			1500 Kg	61KW					
BKQ-1200D-B	1203L	680*1180*1500mm	1450*1990*1710mm	1950 Kg	5KW	380V, 50Hz			
BKQ-1200S-B			2050 Kg	5KW					
BKQ-1500D-B	1500L	680*1180*1870mm	1450*1990*2080mm	2111 Kg	5KW	380V, 50Hz			
BKQ-1500S-B			2211 Kg	5KW					