

## Product Description

**Catalog Number:** 3792

**Product Description:** Costar® 96-well non-treated solid black plate round bottom without lid

**Component Materials:**  
Plate - Virgin Polystyrene, meets *USP, Class VI* requirements for plastic containers and closures. Black concentrate

**Product Dimensions:**

Length of Plate	- 5.035 in.	Diameter of Well @ top	- .270 in.
Width of Plate	- 3.370 in.	Diameter of Well @ bottom	- .250 in.
Depth of Well	- .445 in.	Height without Lid	- .560 in.
Tolerances of Dimensions	- +/- .010 in.	Working Volume per well	- 360 µL

**Surface Characterization:**

Surface is characterized to be medium binding.

**Binding Characteristics:**

The product has been tested for the attribute of well to well protein binding utilizing a fluorescent immunoassay to measure precision. Coefficient of variation must be  $\leq 15\%$ .

**Optical Characteristics:**

The product is made of opaque black polystyrene walls to minimize well to well crosstalk and background fluorescence and /or luminescence.

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes:

Visual examination of the product.

Packaging:

Inspection for accurate labeling, and correct product configuration.

**Lot Number Designation:**

8 Digit Lot Number: First 3 digits - Julian Date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.