

Corning® Gosselin™ Octagonal PET Storage Bottles, Non-sterile

For storage and sampling of solutions

CORNING

Polyethylene terephthalate (PET) bottles are ideal for nonsterile storage and sampling of aqueous solutions. The PET resin makes the bottle lightweight, break-resistant, and as transparent as glass.

Key Features

- ▶ Available in four sizes: 125 mL, 250 mL, 500 mL, and 1L
- ▶ Bulk cases with convenient shrink-wrapped packs of 12 or 24 bottles
- ▶ Break-resistant bottles are ergonomic and an easy alternative to glass
- ▶ Leak-proof, tamper-evident 31.7 mm HDPE screw cap
- ▶ Molded graduations for accurate measurements
- ▶ Bottles are as transparent as glass



Ordering Information

Corning® Gosselin™ Octagonal PET* Storage Bottles, Non-sterile



Molded graduations for accurate measurements



Tamper-evident screw cap

Cat. No.	Description	Qty/Pk	Qty/Cs
P125B-05	125 mL, graduated, non-assembled tamper-evident cap	24	288
P250B-01	250 mL, graduated, non-assembled tamper-evident cap	24	144
P500B-09	500 mL, graduated, non-assembled tamper-evident cap	12	120
P1000B-06	1L, graduated, non-assembled tamper-evident cap	12	48

*Polyethylene terephthalate

Visit www.corning.com/lifesciences for a complete listing of Corning storage bottles and these additional technical resources:

- ▶ Properties of Corning Plasticware
- ▶ Chemical Compatibility of Corning Plasticware

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

CORNING

Corning Incorporated
Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832

China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan

t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6572-9740
f 65 6735-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

CEurope@corning.com

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 020 655 79 28
f 020 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries

t +31 (0) 206 59 60 51
f +31 (0) 206 59 76 73

LATIN AMERICA
grupoLA@corning.com

Brasil
t 55 (11) 3089-7400

Mexico
t (52-81) 8158-8400