



## ROTOR SPECIFICATIONS

<b>Type</b>	Swinging Bucket
<b>Material</b>	Aluminum
<b>Capacity</b>	4 x 2.2 mL
<b>Maximum Speed</b>	55,000 rpm
<b>Tube Dimensions Ø x L</b>	11 x 34 mm
<b>Angle</b>	180°
<b>K Factor</b>	43.8
<b>Net Weight</b>	1.19 kg/2.6 lbs
<b>Warranty</b>	5 years

	<b>RCF</b>	<b>RADIUS</b>
<b>Maximum</b>	258,826 xg	76.6 mm
<b>Average</b>	206,114 xg	61.0 mm
<b>Minimum</b>	153,403 xg	45.4 mm

## ROTOR AND ACCESSORIES ORDERING INFORMATION

## S55-S Rotor Package includes:

<b>Description</b>	<b>Cat. No.</b>
S55-S Rotor	45594
2.2 mL PA Thick-Walled Tubes (qty. 2 pkg. of 50)	45240
Bucket O-ring (qty. 8)	45276
Bucket Stand	45312
Bucket Cap Wrench	45277
Rotor Stand	45275
Lubricant	61556
Vacuum Grease	65937
Instruction Manual	45735

## Accessories:

<b>Description</b>	<b>Cat. No.</b>
Conical Adapter Removal Tool	75000480
Spare Buckets	45358

## LABWARE ORDERING INFORMATION

Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF (xg)	Dims ØxL (mm)	REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
									Type	Cat. No.	Qty. /Pkg.	Description
2.2	2.2	PC Thin-Walled Tube	45315	50	4	55,000	258,826	11 x 34	–	–	–	–
2.2	2.2	PC Ultrathin Tube	75000515	50	4	55,000	258,826	11 x 34	–	–	–	–
2.2	2.2	PA Thick-Walled Tube	45240	50	4	55,000	258,826	11 x 34	–	–	–	–
1.4	1.4	PC Ultracone Tube	75000516	50	4	55,000	258,826	11 x 32	–	75000470	1	Conical Adapter
									–	75000480	1	Conical Adapter Removal Tool
1.4	1.4	PA Ultracone Tube	75000517	50	4	55,000	258,826	11 x 32	–	75000470	1	Conical Adapter
									–	75000480	1	Conical Adapter Removal Tool
1.0	1.0	PC Thick-Walled Tube	45237	100	4	55,000	258,826	11 x 34	–	–	–	–

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

For more information, please visit  
[www.thermoscientific.com/rotors](http://www.thermoscientific.com/rotors).