



The following features are present in the vector based on the nucleotide sequence.

Thymidine kinase (TK) promoter	26-777
<i>Cypridina</i> luciferase gene	790-2451
BGH poly(A) signal.....	2460-2585
SV40 origin/promoter.....	2586-3150
Puromycin resistant gene	3151-3750
SV40 poly(A) signal	3912-3945
beta-lactamase (Amp ^R) gene	4054-4914
pUC replication origin (pUC Ori).....	5093-5897
Transcriptional terminator (Ter)	5898-6505
Lac operator 1 (Lac O1)	6506-6526
Transcriptional pause site (TPS).....	6659-6730