Corning Incorporated Life Sciences

Registered ISO 9001:2008

Product Description

Catalog Number: 9999-132

Product Description: PYREX® Screw Cap, Phenolic, Rubber Liner, G.P.I. 415-13

Component Materials:

Cap - Phenolic (Mineral Fill) Liner - Rubber

Attributes:

Product Dimensions - 16 mm OD x 14 mm High

Autoclaving:

This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes. Caps should be loose on the container to avoid liners being pulled out by the vacuum created as the assembly cools.

Heating / Cooling:

This product is primarily designed for a specific use that does not involve direct heating. In normal use, temperatures will be from -50° to +125°C. This is below the temperature at which thermal failure could be a problem.

Vacuum / Pressure:

This product is a component part and not specifically rated for vacuum or pressure.

Sterilization:

Shipped Not Sterile

Lot Number Designation:

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Revision Date:

16 AUG 2013

Rev No: 1