Corning Incon Life Scier Registere ISO 9001:2	Product Description
Catalog Number:	3335
Product Description:	Costar $\ensuremath{^{ extsf{e}}}$ 6-well cell culture mulitple well plate, with CellBIND $\ensuremath{^{ extsf{e}}}$ surface, flat bottom, with lid
Component Materials: Plate/Lid -	Virgin Polystyrene, meets USP, Class VI requirements for plastic containers and closures
Product Dimensions: Length of Plate Width of Plate Depth of Well Tolerances of Dimensions	- 5.030 in. Diameter of Well @ top - 1.395 in. - 3.365 in. Diameter of Well @ - 1.370 in. - .685 in. Height with Lid - .891 in. - +/050 in. Working Volume per well - 5 mL

Sterilization:

The product is irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137 *Sterilization of healthcare products-Requirements for validation and routine control-Radiation sterilization.* Sterility Assurance Level: SAL 10⁻³

Pyrogens:

The product has been tested and has met the criteria established in the current version of ANSI/AAMI ST 72:20002/(R)2010 *Bacterial Endotoxins - Test methodologies,routine monitoring, and alternative to batch testing.* Results: $\leq 0.1 \text{ EU/mL}$ ($\leq 4 \text{EU/device}$)

Surface Characterization:

Surface is characterized to be hydrophilic and negatively charged, composed of \geq 20% oxygen atoms. This surface composition has been developed for cell attachment and growth.

Cell Attachment and Growth Characteristics:

The product has been tested for the attribute of cell attachment and growth utilizing an attachment- dependent mammalian cell line in a serum supplemented media.

Bovine Spongiform Encephalopathy and Transmissible Spongiform Encephalopathy:

This product is manufactured with animal free materials.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures. Cell Culture Treatment: Visual Attributes: Packaging: Visual and microscopic examination of the product. Inspection for seal and barrier integrity, accurate labeling, and correct product configuration.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 8