

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 8506

Product Description: Corning® 96 Storage Tube Septum Cap, Thermoplastic Elastomer (TPE), Not Sterile, Bulk Packed in Bag, 96 Caps/Bag, 100 Bags/Case (9600 Caps)

Component Materials:

Cap - Thermoplastic Elastomer (PTE)

Attributes:

Cap Dimensions - 8.3 mm OD x 4.2 mm Height

Heating / Cooling:

This product is designed for use at temperatures ranging from -100°C (148°F) to 110°C (230°F).

Sterilization:

Supplied Not Sterile.

Lot Number Designation:

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date; Next 2 digits - Last 2 digits of year of manufacture; Last 3 digits - Batch identification.

Rev No: 1