

**Corning Incorporated
Life Sciences**

Registered
ISO 9001:2008

Product Description

Catalog Number: 1395-45LTC3

Product Description: Corning® Screw Cap, Polypropylene, GL-45, Green

Component Materials:

Cap - Polypropylene (PP), Color: Green.

Attributes:

Product Dimensions - 54 mm OD x 26 mm High

Autoclaving:

This product will withstand repeated autoclave cycles of 121°C at 15 psi for 15 minutes.

Heating / Cooling:

The normal use temperature range for this product is from -50°C to +125°C. The upper temperature limit for the plastic components of this product is 140°C. This product is not designed for direct contact heating.

Vacuum / Pressure:

This product is a component part and not specifically rated for vacuum or pressure.

Sterilization:

Shipped Not Sterile

TSE/BSE:

This product complies with the latest revision of EMEA/410/01, *Note for Guidance on Minimizing the Risk of Transmitting Animal Spongiform Encephalopathy Agents via Human and Veterinary Medicinal Product* by virtue that any component of the product which may include material of bovine origin has been processed per the specific conditions of section 6.4 of EMEA/410/01.

Lot Number Designation:

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Revision Date:

14 FEB 2014

Rev No: 2