



The following features are present in the vector based on the nucleotide sequence.

Multiple cloning site (MCS)	1-42
<i>Cypridina</i> luciferase gene	49-1710
BGH poly(A) signal	1719-1844
SV40 origin/promoter	1845-2409
Puromycin resistant gene	2410-3009
SV40 poly(A) signal	3171-3204
beta-lactamase (Amp ^R) gene	3313-4173
pUC replication origin (pUC Ori)	4352-5156
Transcriptional terminator (Ter)	5157-5764
Lac operator 1 (Lac O1)	5765-5785
Transcriptional pause site (TPS)	5918-5989