

Corning Incorporated  
Life Sciences

Registered  
ISO 9001:2008

## Product Description

**Catalog Number:** 8507

**Product Description:** Corning® 96 Storage Tube Septum Cap, Thermoplastic Elastomer (TPE), Not Sterile, On Mat, 96 Caps/Mat, 10 Mats/Pack, 5 Packs/Case (4800 Caps)

**Component Materials:**

Cap - Thermoplastic Elastomer (PTE)

**Attributes:**

Cap Dimensions - 8.3 mm OD x 4.2 mm Height

Cap Pattern - 8 x 12 on 9.0 mm center

**Heating / Cooling:**

This product is designed for use at temperatures ranging from -100°C (148°F) to 110°C (230°F).

**Sterilization:**

Supplied Not Sterile.

**Lot Number Designation:**

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date; Next 2 digits - Last 2 digits of year of manufacture; Last 3 digits - Batch identification.

Rev No: 2