

Corning® Cryogenic Vials

The advantages are easy to see




CORNING



Corning® Cryogenic Vials

In order to best meet your needs, Corning offers three styles of cryogenic vials as well as storage racks and boxes. Corning also provides a protocol, *Cryogenic Preservation and Storage of Animal Cells* and a technical report, *General Guide for Cryogenically Storing Animal Cell Culture* on our web site at www.corning.com/lifesciences under technical information.

Cryogenic Vials Designs

 <p>External Thread Cryogenic Vial</p> <ul style="list-style-type: none"> ▶ Color-coded polypropylene cap inserts simplify vial identification. Available in variety packs of white, blue, green, red, and yellow. ▶ Silicone washer provides a secure seal. ▶ Easy-to-read black graduations for partial volumes ▶ Self-standing base, self-locking skirt 	 <p>Internal Thread Cryogenic Vial</p> <ul style="list-style-type: none"> ▶ Color-coded polypropylene cap inserts simplify vial identification. Available in variety packs of white, blue, green, red, and yellow. ▶ Silicone washers or rubber O-rings provide a secure seal. ▶ Easy-to-read black graduations for partial volumes ▶ Self-standing base, self-locking skirt 	 <p>External Thread Plug Seal Cap</p> <ul style="list-style-type: none"> ▶ Sure-grip plug seal screw cap ▶ Inner cap ring assures a tight seal.
--	--	---

Warning! Do not use cryogenic vials for storage in the liquid phase of liquid nitrogen. Only store vials in the vapor phase above the liquified gas. Always use appropriate safety equipment when removing vials from cryogenic storage.

Corning External Thread Cryogenic Vials



- ▶ Manufactured from polypropylene to withstand temperatures to -196°C
- ▶ 40% larger marking spot
- ▶ New black graduations
- ▶ Now certified RNase-/DNase-free
- ▶ Vials have a silicone washer for a secure seal
- ▶ Vials may be color coded with inserts (Cat. No. 430499)
- ▶ Self-standing vials have a special base design allowing them to be locked into cryogenic rack and tray (Cat. No. 430525 or 431131) for single-handed manipulation
- ▶ Sterilized by gamma radiation
- ▶ Certified nonpyrogenic

External Thread Cryogenic Vial Ordering Information

Cat. No.	Capacity (mL)	Style	Self-Standing	Qty/Pk	Qty/Cs
430658	1.2	Conical bottom	Yes	50	500
430659	2.0	Round bottom	Yes	50	500
430661	2.0	Round bottom	No	50	500
430662	4.0	Round bottom	Yes	50	500
430663	5.0	Round bottom	Yes	50	500

Corning® Cryogenic Vials



Corning Internal Thread Cryogenic Vials

- ▶ Manufactured from polypropylene to withstand temperatures down to -196°C
- ▶ 40% larger marking spot
- ▶ New black graduations
- ▶ Now certified RNase-/DNase-free
- ▶ Vials have a silicone washer or rubber O-ring for a secure seal
- ▶ Vials may be color coded with inserts (Cat. No. 430499)
- ▶ Self-standing vials have a special base design allowing them to be locked into cryogenic rack and tray (Cat. No. 430525 or 431131) for single-handed manipulation
- ▶ Sterilized by gamma radiation
- ▶ Certified nonpyrogenic

Internal Thread Cryogenic Vial Ordering Information

Cat. No.	Capacity (mL)	Style	Self-Standing	Seal Type	Qty/Pk	Qty/Cs
430487	1.2	Conical bottom	Yes	Washer	50	500
430488	2.0	Round bottom	Yes	Washer	50	500
430489	2.0	Round bottom	No	Washer	50	500
430490	4.0	Round bottom	No	Washer	50	500
430491	4.0	Round bottom	Yes	Washer	50	500
430492	5.0	Round bottom	No	Washer	50	500
430656	5.0	Round bottom	Yes	Washer	50	500
431386	2.0	Round bottom	Yes	Washer	50	250

Corning External Thread Cryogenic Vials with Plug Seal Cap

- ▶ Manufactured from polypropylene to withstand temperatures to -196°C
- ▶ Vials feature an external thread with a traditional plug seal cap design for a secure seal
- ▶ Cap does not accept color-coded inserts
- ▶ Sterilized by gamma radiation
- ▶ Certified nonpyrogenic



External Thread Cryogenic Vial with Plug Seal Cap Ordering Information

Cat. No.	Capacity (mL)	Style	Self-Standing	Qty/Pk	Qty/Cs
430289	2.0	Round bottom	No	50	500

Warning! Do not use cryogenic vials for storage in the liquid phase of liquid nitrogen. Only store vials in the vapor phase above the liquified gas. Always use appropriate safety equipment when removing vials from cryogenic storage.



Corning® Cap Inserts for Cryogenic Vials

- ▶ Caps inserts provide color-coding for easy sample identification.
- ▶ Inserts are packaged in resealable bags.
- ▶ Nonsterile
- ▶ Cap inserts fit all Corning cryogenic vials except Cat. No. 430289.

Cryogenic Cap Inserts Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
430499	Assorted colors, polypropylene cap inserts: 100 each of white, blue, red, green, and yellow	50	500
2015	White polypropylene cap inserts	50	500
2016	Blue polypropylene cap inserts	50	500
2017	Red polypropylene cap inserts	50	500
2018	Green polypropylene cap inserts	50	500
2019	Yellow polypropylene cap inserts	50	500

Corning Cryogenic Vial Racks and Storage Boxes

- ▶ Reusable racks are designed for use with most cryogenic vials.
- ▶ Cat. No. 430525 has a locking feature for use with all Corning self-standing vials.



Cryogenic Vial Rack and Box Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
430525	Polycarbonate rack and tray, holds 30 vials; self-locking design in ice/water bath	1	1
430526	Polycarbonate rack only, holds 30 vials; self-locking design	1	1
431131	Reusable orange polypropylene vial rack, holds 50 vials; self-locking design	2	2
431119	81 Count (9 x 9 array) Cryogenic Box, for 1-2 mL vials	5	10
431120	81 Count (9 x 9 array) Cryogenic Box, for 4-5 mL vials	5	10
431121*	100 Count (10 x 10 array) Cryogenic Box, for 1-2 mL vials	5	10

*431121 accepts internally threaded cryogenic vials only.

For additional product, accessory, and application information, please visit www.corning.com/lifesciences, or call 1.800.492.1110. International customers, please call 1.978.442.2200 or contact your local support office (listed below).

CORNING

Corning Incorporated Life Sciences

Tower 2, 4th Floor
900 Chelmsford St.
Lowell, MA 01851
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

Worldwide Support Offices

ASIA / PACIFIC

Australia/New Zealand
t 0402-794-347

China
t 86 21 2215 2888
f 86 21 6215 2988

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6733-6511
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries

t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

LATIN AMERICA

Brasil
t (55-11) 3089-7419
f (55-11) 3167-0700

Mexico
t (52-81) 8158-8400
f (52-81) 8313-8589