

**Corning Incorporated  
Life Sciences**

Registered  
ISO 9001:2008

## Product Description

**Catalog Number:** 3968

**Product Description:** Corning® CellSTACK™ - 33mm orange vent cap, 0.2um membrane

**Component Materials:**

- Cap - Virgin High Density Polyethylene, meets *USP, Class VI* requirements for plastic containers and closures. Heavy metal free (meets CONEG requirements) color concentrate.
- Filter - 0.2µm polyester non-woven membrane, meets *USP, Class VI* requirements for plastic containers and closures.

**Product Dimensions:**

|            |   |              |        |   |          |
|------------|---|--------------|--------|---|----------|
| Diameter   | - | 1.460 in     | Height | - | .840 in. |
| Tolerances | - | +/- .050 in. |        |   |          |

**Sterilization:**

The lot has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137 (AAMI TIR 27) *Sterilization of health care products-Requirements for validation and routine control-Radiation sterilization*

Sterility Assurance Level: SAL 10<sup>-6</sup>

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes: Visual examination of the product.

Packaging: Inspection for seal and barrier integrity, accurate labeling, and correct product configuration.

**Lot Number Designation:**

8 Digit Lot Number: First 3 digits - Julian date, start of manufacture; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 1