

**Corning Incorporated
Life Sciences**

**Registered
ISO 9001:2008**

Product Description

Catalog Number: 6509

Product Description: Costar® 96-well Thermowell™ plate Model P

Component Materials:

Plate - Polycarbonate

Product Dimensions:

Length of Plate	-	4.62 in.	Diameter of Well @ top	-	.254 in.
Width of Plate	-	3.12 in.	Depth of Well	-	.520 in.
Tolerances of Dimensions	-	+/- .050 in.	Total Volume per well	-	200µL

Functional Characteristics:

The product is thermally stable up to 135° C.

Compatible Thermal Cyclers:

Perkin Elmer GeneAmp PCR System 9600, Barnstead Thermolyne Amplitron II, Techne Cyclogene and Gene E with 96 x 0.2 mL block.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes: Visual examination of the product.

Packaging: Inspection accurate labeling, and correct product configuration.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 2