

**Corning Incorporated
Life Sciences**

**Registered
ISO 9001:2008**

Product Description

Catalog Number: 6542

Product Description: Thermowell ® 0.2ml PCR tube strips 1X8, clear

Component Materials:
Tube - Polypropylene

Product Dimensions:
Overall tube length w/out cap - .815" I.D. opening - .211"
Tolerances - +/- .005"

RNase/DNase Testing:
This product has been tested and is free of any detectable RNase/DNase contamination.

Performance Testing:
Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.
Visual Attributes: Visual examination of the product.
Packaging: Inspection for accurate labeling, and correct product configuration.

Lot Number Designation:
8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 1