

Corning Incorporated  
Life Sciences

Registered  
ISO 9001:2008

## Product Description

**Catalog Number:** 4515-36L

**Product Description:** Direct Drive Paddle Assembly, 36 L

**Component Materials:**

Blade, Holder, Sleeve, Plug - Polytetrafluoroethylene (PTFE)  
Coupler, Screws - Stainless Steel

**Attributes:**

Product Dimensions - 146 mm Wide x 433 mm High without Coupler

**Autoclaving:**

This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes.

**Heating / Cooling:**

This product is primarily designed for a specific use that does not involve direct heating. In normal use, temperatures will be from -50° to +125°C. This is below the temperature at which thermal failure could be a problem.

**Vacuum / Pressure:**

This product is not designed for and is not rated for vacuum or pressure applications.

**Sterilization:**

Shipped Not Sterile

**Lot Number Designation:**

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

**Revision Date:**

5 AUG 2013

Rev No: 1