

Corning Incorporated  
Life Sciences

Registered  
ISO 9001:2008

# Product Description

**Catalog Number:** 431413

**Product Description:** Glass Microfiber Prefilters, 79 mm Square, for Use with 1000 mL Filter Funnels

**Component Materials:**  
0.32 mm Thick Binder-Free Microfiber Glass, Grade 161

**Attributes:**  
Product Dimensions - 78.0 mm Square x 0.32 mm Thick

**Lot Number Designation:**  
Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 1