ViralSEQ[™] Mouse Minute Virus (MMV) Real-Time PCR Detection Kit

Catalog Number 4444415

Pub. No. 4444406 Rev. G

Note: For safety and biohazard guidelines, see the "Safety" appendix in the *ViralSEQ* Mouse Minute Virus (MMV) Real-Time PCR Detection Kit User Guide (Pub. No. 4445235). Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Product description

The ViralSEQ $^{\text{\tiny{M}}}$ Mouse Minute Virus (MMV) Real-Time PCR Detection Kit is designed for routine research screening and detection of MMV. This assay, which is based on TaqMan $^{\text{\tiny{M}}}$ technology, has been highly optimized for sensitive and specific detection of MMV and, when coupled with the recommended sample preparation, provides the sensitivity needed for routine in-process testing during product manufacturing.

For more information about reagent use, refer to the $ViralSEQ^{TM}$ Mouse Minute Virus (MMV) Real-Time PCR Detection Kit User Guide (Pub. No. 4445235).

Contents and storage

| Contents | Cap color | Amount | Storage | | |
|---|-----------|--------------|---------------------------------------|--|--|
| Box 1, SEQ Real-Time PCR Core Kit | | | | | |
| 2X Environmental Master Mix 2.0 | | 2 × 0.825 mL | Protect from light. | | |
| Negative Control (water) | | 1.0 mL | −25°C to −15°C | | |
| Box 2, MMV Real-Time PCR Assay Mix | | | | | |
| 10× MMV Assay Mix | | 0.300 mL | Protect from light. -25°C to -15°C | | |
| Box 3, MMV Discriminatory Positive/Extraction Control | | | | | |
| MMV Positive Control, 1,000 copies/μL | | 0.700 mL | −25°C to −15°C | | |

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.





Life Technologies Ltd | 7 Kingsland Grange | Woolston, Warrington WA1 4SR | United Kingdom For descriptions of symbols on product labels or product documents, go to **thermofisher.com/symbols-definition**.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. 4444406 (English)

| Revision | Date | Description | |
|----------|-----------------|--|--|
| G | 4 February 2019 | Update template and legal information. | |
| F | November 2012 | Baseline for this revision history. | |

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license.

