

Superspeed Rotors

Thermo Scientific Fiberlite F12-6x500 LEX Rotor



ROTOR SPECIFICATIONS

Type	Fixed Angle	
Material	Carbon Fiber Composite	
Capacity	6 x 500 mL	
Maximum Speed	12,000 rpm	
Tube Dimensions Ø x L	70 x 160 mm	
Angle	20°	
K Factor	2,302	
Net Weight	9.1 kg/20 lbs	
Warranty	15 years	
Certifications		
	RCF	RADIUS
Maximum	24,471 x g	152 mm
Minimum	6,601 x g	41 mm

ROTOR AND ACCESSORIES ORDERING INFORMATION

		THERMO SCIENTIFIC	
		SORVALL LYNX 6000	SORVALL LYNX 4000
Fiberlite F12-6x500 LEX	Max. Speed	12,000 rpm	12,000 rpm
	Max. RCF	24,471 x g	24,471 x g

Fiberlite F12-6x500 LEX Rotor Package includes:

Description	Cat. No.
Fiberlite F12-6x500 LEX Rotor	096-062375
500 mL Fiberlite PPCO Bottle with PPGF Cap and PP Plug (qty. 6)	010-1493
Rotor Lid Assembly	099-062375
O-Ring Refresher Kit ¹	021-062375

LABWARE ORDERING INFORMATION

Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF ² (x g)	Dims ØxL (mm)	REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
									Type	Cat. No.	Qty. /Pkg.	Description
500.0	500.0	Fiberlite PPCO Oak Ridge Bottle	010-1493	6	6	12,000	24,500	70 x 160	Closure	Included	6	PPGF Cap w/ PP Plug
500.0	500.0	Fiberlite PC Oak Ridge Bottle	010-1494	6	6	12,000	24,500	70 x 160	Closure	Included	6	PPGF Cap w/ PP Plug
500.0	400.0	Nalgene PPCO Bottle	3141-0500	24	6	9,000	13,700	70 x 160	Closure	Included	24	PP Sealing
500.0	400.0	Nalgene PC Bottle	3140-0500	24	6	9,000	13,700	70 x 160	Closure	Included	24	PP Sealing
250.0	250.0	Fiberlite PPCO Oak Ridge Bottle	010-1495	6	6	12,000	22,260	61 x 124	Closure	Included	6	PPGF Cap w/ PP Plug
									Adapter	010-0151	2	1 place/adaptor

Biocontainment certification by Public Health England, Porton Down, UK.

¹ Includes 5 O-Rings and 3 O-Ring Grease Packets.

Refer to manufacturer's recommendations for product performance information.

LABWARE ORDERING INFORMATION CONTINUED

Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF ² (x g)	Dims ØxL (mm)	REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
									Type	Cat. No.	Qty. /Pkg.	Description
250.0	250.0	Fiberlite PC Oak Ridge Bottle	010-1496	6	6	12,000	22,260	61 x 124	Closure	Included	6	PPGF Cap w/ PP Plug
									Adapter	010-0151	2	1 place/adapter
250.0	250.0	Nalgene PPCO Bottle	3141-0250	36	6	12,000	22,260	61 x 133	Closure	Included	36	PP Sealing
									Adapter	010-0151	2	1 place/adapter
250.0	250.0	Nalgene PC Bottle	3140-0250	36	6	12,000	22,260	61 x 135	Closure	Included	36	PP Sealing
									Adapter	010-0151	2	1 place/adapter
250.0	–	Corning™ Disposable Conical Bottle	–	–	6	–	–	–	Adapter	010-1135	2	1 place/adapter
175.0	–	Nalgene PPCO Conical Wide-Mouth Bottle	3143-0175	36	6	–	–	61 x 144	Closure	Included	36	PP Sealing
									Adapter	010-0152	2	1 place/adapter
175.0	–	Nalgene PC Conical Wide-Mouth Bottle	3144-0175	36	6	–	–	61 x 144	Closure	Included	36	PP Sealing
									Adapter	010-0152	2	1 place/adapter
80.0	73.0	PC Oak Ridge Tube	010-0515	6	6	12,000	19,190	38 x 109	Closure	Included	6	PP Sealing
									Adapter	010-1114	2	1 place/adapter
80.0	76.0	PPCO Oak Ridge Tube	010-1280	6	6	12,000	19,190	38 x 109	Closure	Included	25	PP Screw Top
									Adapter	010-1114	2	1 place/adapter
50.0	–	Nunc™ Disposable Conical Tube	339653	25	6	–	–	30 x 115	Closure	Included	25	PP Sealing
									Adapter	010-1102	2	1 place/adapter
50.0	–	Corning Disposable Conical Tube	–	–	6	–	–	–	Adapter	010-1102	2	1 place/adapter
50.0	–	Falcon™ Disposable Conical Tube	–	–	6	–	–	–	Adapter	010-1102	2	1 place/adapter
50.0	46.0	Nalgene PC Oak Ridge Tube	3138-0050	50	12	12,000	22,580	29 x 108	Closure	Included	50	PP Sealing
									Adapter	010-1112	2	2 places/adapter
50.0	46.0	Nalgene PPCO Oak Ridge Tube	3139-0050	50	12	12,000	22,580	29 x 108	Closure	Included	50	PP Sealing
									Adapter	010-1112	2	2 places/adapter
50.0	46.0	Nalgene PSF Oak Ridge Tube	3137-0050	50	12	12,000	22,580	29 x 108	Closure	Included	50	PP Sealing
									Adapter	010-1112	2	2 places/adapter
50.0	46.0	Nalgene FEP Oak Ridge Tube	3114-0050	10	12	12,000	22,580	29 x 108	Closure	Included	10	ETFE Screw Cap
									Adapter	010-1112	2	2 places/adapter
30.0	30.0	Nalgene PC Oak Ridge Tube	3138-0030	50	18	12,000	20,640	26 x 102	Closure	Included	50	PP Sealing
									Adapter	010-1115	2	3 places/adapter
30.0	30.0	Nalgene PPCO Oak Ridge Tube	3139-0030	50	18	12,000	20,640	26 x 102	Closure	Included	50	PP Sealing
									Adapter	010-1115	2	3 places/adapter
30.0	30.0	Nalgene PSF Oak Ridge Tube	3137-0030	50	18	12,000	20,640	26 x 102	Closure	Included	50	PP Sealing
									Adapter	010-1115	2	3 places/adapter
30.0	30.0	Nalgene FEP Oak Ridge Tube	3114-0030	10	18	12,000	20,640	26 x 102	Closure	Included	10	ETFE Screw Cap
									Adapter	010-1115	2	3 places/adapter
16.0	16.0	Nalgene PC Oak Ridge Tube	3138-0016	50	36	12,000	20,640	18 x 107	Closure	Included	50	PP Sealing
									Adapter	010-1105	2	14 places/adapter
16.0	16.0	Nalgene PPCO Oak Ridge Tube	3139-0016	50	36	12,000	20,640	18 x 107	Closure	Included	50	PP Sealing
									Adapter	010-1105	2	14 places/adapter
15.0	15.0	Corning Disposable Conical Tube	–	–	36	–	–	–	Adapter	010-1099	2	12 places/adapter
15.0	15.0	Falcon Disposable Conical Tube	–	–	36	–	–	–	Adapter	010-1099	2	12 places/adapter
10.0	8.0	Nalgene PPCO Oak Ridge Tube	3139-0010	50	42	12,000	20,640	16 x 100	Closure	Included	50	PP Sealing
									Adapter	010-1308	2	7 places/adapter
10.0	8.0	Nalgene PC Oak Ridge Tube	3138-0010	50	42	12,000	20,640	16 x 100	Closure	Included	50	PP Sealing
									Adapter	010-1308	2	7 places/adapter
10.0	–	BD Vacutainer™ Tube	–	–	42	–	–	16 x 100	Adapter	010-1103	2	7 places/adapter
3.0	–	BD Vacutainer Tube	–	–	84	–	–	10 x 64	Adapter	010-1137	2	14 places/adapter

Refer to manufacturer's recommendations for product performance information.

Find out more at thermofisher.com/rotors

© 2012, 2013, 2014, 2021 Thermo Fisher Scientific Inc. All rights reserved. Falcon and Corning are registered trademarks of Corning. BD Vacutainer is a registered trademark of Becton, Dickinson and Company. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

ThermoFisher
SCIENTIFIC