CORNING

Ham's F-10 Medium

Formulation

Cat. No.	10-070
Description	Liquid, 1x
Units	mg/L
Components	
Inorganic Salts	
CaCl ₂ (anhydrous)	33.30
CuSO ₄ (anhydrous)	0.0016
FeSO ₄ • 7H ₂ O	0.834
КСІ	285.00
KH ₂ PO ₄	83.00
MgSO ₄ (anhydrous)	74.60
NaCl	7400.00
Na ₂ HPO ₄ (anhydrous)	153.70
NaHCO ₃	1200.00
ZnSO ₄ • 7H ₂ O	0.0288
Amino Acids	
L-Alanine	9.00
L-Arginine • HCl	211.00
L-Asparagine • H_2O	15.00
L-Aspartic acid	13.00
L-Cysteine • HCl • H ₂ O	35.10
L-Glutamic acid	14.70
L-Glutamine	146.20
Glycine	7.50
L-Histidine • HCl • H_2O	23.00
L-Isoleucine	2.60
L-Leucine	13.00
L-Lysine • HCl	29.00
L-Methionine	4.48
L-Phenylalanine	5.00
L-Proline	11.50
L-Serine	10.50
L-Threonine	3.57
L-Tryptophan	0.60
L-Tyrosine • 2Na • 2H ₂ O	2.61
L-Valine	3.50
Vitamins	
Biotin	0.024
D-Calcium pantothenate	0.715
Choline chloride	0.698
Folic acid	1.32
<i>i</i> -Inositol	0.541
Nicotinamide	0.615
Pyridoxine • HCl	0.206
Riboflavin	0.376
Thiamine • HCl	1.00
Vitamin B ₁₂	1.36

Components (continued)	(10-070)
Other	
D-Glucose	1100.00
Hypoxanthine, Na	4.77
DL-Thioctic (lipoic) acid	0.20
Phenol red • Na	1.20
Sodium pyruvate	110.00
Thymidine	0.73



At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

CORNING | FALCON' CELIGIO AXYGEN, PYREX' GOSSELIN'

For additional product or technical information, please visit **www.cellgro.com** or call 1.800.235.5476.

Mediatech, Inc. *A Corning Subsidiary* 9345 Discovery Boulevard Manassas, VA 20109 t 800.235.5476 t 703.471.5955 f 703.467.9851

www.cellgro.com

Corning cellgro products are manufactured by Mediatech, Inc. For a listing of Corning trademarks, visit us at www.cellgro.com/about-us/trademarks.html All other trademarks in this document are the property of their respective owners.