

Corning Incorporated  
Life Sciences

Registered  
ISO 9001:2008

# Product Description

**Catalog Number:** 1395-25SS

**Product Description:** Corning® Septa, Silicone, for GL25 Open Top PBT Screw Cap

**Component Materials:**  
Silicone Rubber

**Attributes:**  
Product Dimensions - 23.1 mm OD x 2.5 mm Thick

**Autoclaving:**  
This product will withstand repeated autoclave cycles of 121°C at 15 psi for 15 minutes.

**Heating / Cooling:**  
The normal use temperature range for this product is from -50°C to +125°C. This product is not designed for direct contact heating.

**Vacuum / Pressure:**  
This product is a component part and not specifically rated for vacuum or pressure.

**Sterilization:**  
Shipped Not Sterile

**TSE/BSE:**  
Materials of animal origin are not used in the manufacture of this product.

**Lot Number Designation:**  
Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

**Revision Date:**  
11 SEPT 2014

Rev No: 1