

Hot Air Sterilizer

Application:

Hot air sterilizer is used to sterilize small surgical instrument, glass, petri dishes etc.

This device is suitable for medical and veterinary clinics, hospitals, pharmacies, health care centers, and laboratories.



HAS-T25/HAS-T50/HAS-T105/HAS-T200

HAS-T70

Features:

- * Adjustable sterilization time.
- * Top air outlet with adjustable air volume.
- * Over-temperature alarm.
- * Heating up fast, forced-air convection, fast sterilization.
- * Fan inside, fully sealed inside, no air circulation with external air, temperature rise fast.
- * High precision microprocessor temperature control; 304 stainless steel inner chamber.

Technical Parameters:

Model	HAS-T25	HAS-T50	HAS-T70	HAS-T105	HAS-T200
Capacity	25L	50L	70L	105L	200L
Temp. Range	RT+10-250°C				
Temp. Precision	0.1°C				
Temp. Fluctuation	±0.5°C				
Ambient Temp.	5-40°C				
Timing Range	1-9999min				
Shelves No.	2 pcs				
Consumption	500W	750W	1050W	1500W	2000W
Power Supply	AC110/220V±10%, 50/60Hz				
Internal Size(W*D*H) mm	300*300*270	420*350*350	400*375*470	550*350*520	800*550*570
External Size(W*D*H) mm	590*470*440	710*520*520	540*550*870	840*540*720	890*730*770
Package Size(W*D*H) mm	710*600*550	810*840*800	670*660*1040	950*680*880	1020*890*930
Gross Weight(kg)	45	50	75	80	100