

TECHNICAL SHEET

PRODUCT NAME: Anti-Influenza A mAb (clone Y77)

Code: MT-18Y77

Clone: Y77

Isotype: IgG₁

Physical State: Liquid

Host Animal (source): Mouse

Specificity: Y77 is specific for Influenza A detection.

Cross Reactivity: No cross reactivity with *Adenovirus, Respiratory Syncytial Virus.*

Purification: Anti-Influenza A mAb (clone Y77) has been purified in a multi-

step procedure.

≥95% pure by SDS-PAGE.

Storage Conditions: Storage at 2-8°C.

Health & Safety Information: This product contains sodium azide as preservative. The amount

is 0.09%, this substance is poisonous and hazardous and it should

be handled by trained staff only.

Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance

with the respective national regulations.

Presentation: Citrate pH: 6, 0.09% Sodium Azide (NaN₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 18/01/2019



www.certest.es