

A Complete Pipetting Solution

Axygen® Axypet® Pro Single- and Multi-channel Pipettors and Axxygen MultiRack Pipet Tip System



Axygen Axypet Pro pipettors and the Axxygen MultiRack tip system have been engineered to provide the highest levels of comfort, accuracy, and precision.

The lightweight construction, contoured handle, and 4-digit counter of Axxygen Axypet Pro pipettors ensure comfortable pipetting. All feature smooth plunger movement and extremely low pipetting forces to reduce wrist strain and fatigue. Pushbuttons on the pipettors are color-coded to match the racks of Axxygen MultiRack tips to enable easy size identification and reduce the risk of cross-contamination.

The Axxygen MultiRack tip system is an interchangeable, complete system offered in bulk, racked, filter, and reload packaging styles to meet the variety of needs in the lab setting. The MultiRack tip system offers standard and extended-length tips allowing volumes from 0.1 μL to 1250 μL . For racked and reload packaging styles, the tips are pre-loaded into easy-load, color-coded deck inserts. With a small benchtop footprint and reusable rack, this system minimizes environmental impact and lab space. The stackable design saves critical bench space. Axxygen MultiRack tip system products are made in the USA and meet ISO standards. All tips are free of detectable RNase, DNase, DNA, and pyrogens.

Axygen Axypet Pro Pipettors

- ▶ Ergonomic handle design fits both the right and left hand
- ▶ Convenient one-hand volume setting with auto-lock system prevents accidental volume change during work
- ▶ 4-digit counter provides greater precision and is positioned for visibility at all times
- ▶ Single-channel pipettors are compatible with very narrow tubes
- ▶ Full autoclavability and UV-resistance helps prevent sample contamination
- ▶ The multi-channel shafts retract individually for easier tip loading and ejection
- ▶ Unique ejector lever system enables effortless tip ejection
- ▶ Easy in-lab calibration helps you quickly adjust your pipettor to various liquid types
- ▶ Individually tested and supplied with a certificate of quality

Axygen MultiRack Pipet Tips

- ▶ Reusable rack
- ▶ Small footprint and stackable design for maximum bench space
- ▶ Extended tips offered in 10, 200, 300, and 1000 μL sizes (1000 μL extended length allows up to 1250 μL volume)

Compatibility Chart



Axygen Axypet Pro Pipettor		Axygen MultiRack Pipet Tips									
		Non-filtered				Filtered					
Cat. No.	Description	10 µL	200 µL	300 µL	1000 µL	10 µL	100 µL	200 µL	300 µL	1000 µL	
AP-2-P	0.1-2 µL single-channel pipettor	■				■					
AP-10-P	0.5-10 µL single-channel pipettor	■				■					
AP-20-P	2-20 µL single-channel pipettor		■								
AP-50-P	5-50 µL single-channel pipettor		■				■				
AP-100-P	10-100 µL single-channel pipettor		■				■				
AP-200-P	20-200 µL single-channel pipettor		■					■			
AP-1000-P	100-1,000 µL single-channel pipettor				■					■	
AP-8-10-P	0.5-10 µL, 8-channel pipettor	■				■					
AP-8-50-P	5-50 µL, 8-channel pipettor		■				■				
AP-8-200-P	10-200 µL, 8-channel pipettor			■					■		
AP-8-300-P	30-300 µL, 8-channel pipettor			■					■		
AP-12-10-P	0.5-10 µL, 12-channel pipettor	■				■					
AP-12-50-P	5-50 µL, 12-channel pipettor		■				■				
AP-12-200-P	10-200 µL, 12-channel pipettor			■					■		
AP-12-300-P	30-300 µL, 12-channel pipettor			■					■		

Ordering Information

Axygen® Axypet® Pro Pipettors

Single-channel

Cat. No.	Volume Range (µL)	Accuracy (%)	Precision (%)	Non-filtered Tip (µL)	Qty/Cs
AP-2-P	0.1 - 2	±40.0 to ±1.5	≤12.0 to ≤0.7	10	1
AP-10-P	0.5 - 10	±4.0 to ±1.0	≤2.8 to ≤0.4	10	1
AP-20-P	2 - 20	±3.5 to ±0.8	≤1.5 to ≤0.3	200	1
AP-50-P	5 - 50	±3.5 to ±0.8	≤1.3 to ≤0.3	200	1
AP-100-P	10 - 100	±3.0 to ±0.8	≤1.0 to ≤0.2	200	1
AP-200-P	20 - 200	±2.0 to ±0.6	≤0.7 to ≤0.2	200	1
AP-1000-P	100 - 1000	±2.5 to ±0.6	≤0.6 to ≤0.2	1000	1

8-channel

AP-8-10-P	0.5 - 10	±10.0 to ±2.0	≤8.0 to ≤1.2	10	1
AP-8-50-P	5 - 50	±4.0 to ±1.6	≤2.5 to ≤0.6	200	1
AP-8-200-P	10 - 200	±5.0 to ±1.0	≤3.0 to ≤0.6	200/300	1
AP-8-300-P	30 - 300	±3.0 to ±1.0	≤1.5 to ≤0.6	300	1

12-channel

AP-12-10-P	0.5 - 10	±10.0 to ±2.0	≤8.0 to ≤1.2	10	1
AP-12-50-P	5 - 50	±4.0 to ±1.6	≤2.5 to ≤0.6	200	1
AP-12-200-P	10 - 200	±5.0 to ±1.0	≤3.0 to ≤0.6	200/300	1
AP-12-300-P	30 - 300	±3.0 to ±1.0	≤1.5 to ≤0.6	300	1

Axygen Axypet Pro Starter Kit

Cat. No.	Description	Qty/Cs
AP-STR-KIT-P	Axygen Axypet Pro Starter Kit Includes: 4 single-channel pipettors: 0.5 to 10 µL (Cat. No. AP-10-P), 2 to 20 µL (Cat. No. AP-20-P), 20 to 200 µL (Cat. No. AP-200-P), 100 to 1000 µL (Cat. No. AP-1000-P); universal linear stand for four single-channel pipettors (Cat. No. AP-USD4-P); Axygen MultiRack pipet tips: 10 µL (Cat. No. MR-10XT-R), 200 µL (Cat. No. MR-200-R), 1000 µL (Cat. No. MR-1000XT-R)	1

Accessories

AP-USD	Pipettor stand for one single-channel, 8-channel, or 12-channel pipettor	1
AP-USD4-P	Universal linear stand for four pipettors, red (four single-channel, maximum of two 8-channel or 12-channel pipettors)	1
AP-USD4	Universal linear stand for four pipettors, transparent (four single-channel, maximum of two 8-channel or 12-channel pipettors)	1
AP-USD6	Universal linear stand for six pipettors, transparent (six single-channel, maximum of three 8-channel or 12-channel pipettors)	1



Ordering Information

Axygen® MultiRack Pipet Tips

Bulk

Cat. No.	Description	Quantity
MR-10XT	10 µL MultiRack extended length pipet tip, bulk	1,000/pk, 10 pks/cs
MR-200	200 µL MultiRack pipet tip, bulk	1,000/pk, 10 pks/cs
MR-1000XT	1000 µL MultiRack extended length pipet tip, bulk	768/pk, 10 pks/cs

Racked

MR-10XT-R	10 µL MultiRack extended length pipet tip, racked	960/pk, 5 pks/cs
MR-200-R	200 µL MultiRack pipet tip, racked	960/pk, 5 pks/cs
MR-300-R	300 µL MultiRack pipet tip, racked	960/pk, 5 pks/cs
MR-300XT-R	300 µL MultiRack pipet tip, extended length, racked	768/pk, 5 pks/cs
MR-1000XT-R	1000 µL MultiRack extended length pipet tip, racked	768/pk, 5 pks/cs

Racked, Sterile

MR-10XT-R-S	10 µL MultiRack extended length pipet tip, racked, sterile	960/pk, 5 pks/cs
MR-200-R-S	200 µL MultiRack pipet tip, racked, sterile	960/pk, 5 pks/cs
MR-300-R-S	300 µL MultiRack pipet tip, racked, sterile	960/pk, 5 pks/cs
MR-300XT-R-S	300 µL MultiRack extended length pipet tip, racked, sterile	768/pk, 5 pks/cs
MR-1000XT-R-S	1000 µL MultiRack extended length pipet tip, racked, sterile	768/pk, 5 pks/cs

Axygen MultiRack Reload System

Cat. No.	Description	Quantity
MR-10XT-RL	10 µL MultiRack extended length pipet tip, reload	960/pk, 4 pks/cs
MR-200-RL	200 µL MultiRack pipet tip, reload	960/pk, 4 pks/cs
MR-300-RL	300 µL MultiRack pipet tip, reload	960/pk, 4 pks/cs
MR-1000XT-RL	1000 µL MultiRack extended length pipet tip, reload	960/pk, 4 pks/cs

Axygen MultiRack Filter Pipet Tips


Racked, Filter, Sterile

Cat. No.	Description	Quantity
MRF-10XT-R-S	10 µL MultiRack extended length filtered pipet tip, racked, sterile	960/pk, 5 pks/cs
MRF-20-R-S	20 µL MultiRack filtered pipet tip, racked, sterile	960/pk, 5 pks/cs
MRF-100-R-S	100 µL MultiRack filtered pipet tip, racked, sterile	960/pk, 5 pks/cs
MRF-200NX-R-S	200 µL NX Axxygen MultiRack filtered pipet tip, racked, sterile	960/pk, 5 pks/cs
MRF-200XT-R-S	200 µL MultiRack extended length filtered pipet tip, racked, sterile	768/pk, 5 pks/cs
MRF-300XT-R-S	300 µL MultiRack extended length filtered pipet tip, racked, sterile	768/pk, 5 pks/cs
MRF-1000XT-R-S	1000 µL MultiRack extended length, filtered pipet tip, racked, sterile	768/pk, 5 pks/cs

Axygen MultiRack Maxymum Recovery® Pipet Tips

Racked, Sterile, Axxygen Maxymum Recovery Surface

Cat. No.	Description	Quantity
MR-10XT-L-R-S	10 µL MultiRack extended length pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MR-200-L-R-S	200 µL MultiRack natural pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MR-300-L-R-S	300 µL MultiRack pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MR-300XT-L-R-S	300 µL MultiRack extended length pipet tip, Maxymum Recovery surface, racked, sterile	768/pk, 5 pks/cs
MR-1000XT-L-R-S	1000 µL MultiRack extended length pipet tip, Maxymum Recovery surface, racked, sterile	768/pk, 5 pks/cs



Ordering Information (Continued)

Racked, Filter, Sterile, Axygen® Maxymum® Recovery Surface

Cat. No.	Description	Quantity
MRF-10XT-L-R-S	10 µL MultiRack extended length filtered pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MRF-20-L-R-S	20 µL MultiRack filtered pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MRF-100-L-R-S	100 µL MultiRack filtered pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MRF-200NX-L-R-S	200 µL NX MultiRack filtered pipet tip, Maxymum Recovery surface, racked, sterile	960/pk, 5 pks/cs
MRF-200XT-L-R-S	200 µL MultiRack extended length filtered pipet tip, Maxymum Recovery surface, racked, sterile	768/pk, 5 pks/cs
MRF-300XT-L-R-S	300 µL MultiRack extended length filtered pipet tip, Maxymum Recovery surface, racked, sterile	768/pk, 5 pks/cs
MRF-1000XT-L-R-S	1000 µL MultiRack extended length, filtered pipet tip, Maxymum Recovery surface, racked, sterile	768/pk, 5 pks/cs

Axygen MultiRack Tip Station

Cat. No.	Description	Quantity
MR-10XT-TS	10 µL MultiRack extended length pipet tip station	192/layer, 5 layers/stack, 960/stack, 5 stacks/cs
MR-10XT-TS-S	10 µL MultiRack extended length pipet tip station, sterile	192/layer, 5 layers/stack, 960/stack, 5 stacks/cs
MR-200-TS	200 µL MultiRack pipet tip station	192/layer, 5 layers/stack, 960/stack, 5 stacks/cs
MR-200-TS-S	200 µL MultiRack pipet tip station, sterile	192/layer, 5 layers/stack, 960/stack, 5 stacks/cs

Accessories

Cat. No.	Description	Quantity
MR-R-A	Empty rack for Axygen MultiRack 10 µL extended length, 200 µL, 200 µL NX, 300 µL	5/pk, 1 pk/cs
MR-R-B	Empty rack for Axygen MultiRack 200 µL extended length, 300 µL extended length, and 1000 µL extended length	5/pk, 1 pk/cs



For more specific information on claims, visit www.corning.com/certificates.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

CORNING

Corning Incorporated
Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

ASIA/PACIFIC
Australia/New Zealand
t 61 427286832

Chinese Mainland
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6572-9740
f 65 6735-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE
CEurope@corning.com

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 020 655 79 28
f 020 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries
t +31 (0) 206 59 60 51
f +31 (0) 206 59 76 73

LATIN AMERICA
grupoLA@corning.com

Brazil
t 55 (11) 3089-7400

Mexico
t (52-81) 8158-8400