

Ordering Information

Cat#	CKS-L010Q125	CKS-L015Q125	CKS-L020Q125	CKS-L025Q125	CKS-L030Q125	CKS-L035Q125	
Unit parameter							
Geometric capacity (L)	10	15	20	25	30	35	
Neck diameter(mm)	125	125	125	125	125	125	
Static evaporation loss rate (L/Day)	0.42	0.36	0.37	0.37	0.33	0.36	
Static holding time (Day)	24	42	54	67	90	97	
Storage capacity							
Racks and vial boxes	Number of racks (ea)	1	1	1	1	6/7	6/7
	Diameter of racks (mm)	82x84	82x84	82x84	82x84	82x84	82x84
	Diameter of boxes (mm)	76x76	76x76	76x76	76x76	76x76	76x76
	Number of layers per rack (ea)	4	3	4	5	4	5
	Number of vials (ea)	100	75	100	125	600/700	750/875
25ml (single)	Number of racks (ea)	—	—	—	—	18	18
	Number of layers per rack (ea)	—	—	—	—	2	2
	Number of blood bags (ea)	—	—	—	—	36	36
25ml (multi)	Number of racks (ea)	—	—	—	—	7	7
	Number of layers per rack (ea)	—	—	—	—	2	2
	Number of boxes per						

25ml (multi)	Number of boxes per layer (ea)	—	—	—	—	4	4
	Number of blood bags (ea)	—	—	—	—	56	56
50ml (single)	Number of racks (ea)	—	—	—	—	—	—
	Number of layers per rack (ea)	—	—	—	—	—	—
	Number of blood bags (ea)	—	—	—	—	—	—
50ml (multi)	Number of racks (ea)	—	—	—	—	—	—
	Number of layers per rack (ea)	—	—	—	—	—	—
	Number of boxes per layer (ea)	—	—	—	—	—	—
	Number of blood bags (ea)	—	—	—	—	—	—
Containers							
Height (mm)	554	589	670	700	705	748	
Outer diameter (mm)	300	394	394	394	461	461	
Empty weight (kg)	5.9	8.5	9.7	10.2	12.9	14.2	
Accessories included							
Lockable lid	√	√	√	√	√	√	
Protective holding bag	√	√	√	√	√	√	
Optional accessories							

Optional accessories						
Intelligent monitor	√	√	√	√	√	√
Carrying cart	—	—	—	—	√	√

Col#	CKS-L047Q127	CKS-L050Q125	CKS-L065Q216	CKS-L095Q216	CKS-L115Q216	CKS-L145Q216
Unit parameter						
Geometric capacity (L)	47	50	65	95	115	145
Neck diameter(mm)	125	125	216	216	216	216
Static evaporation loss rate (L/Day)	0.36	0.42	0.78	0.97	0.94	0.96
Static holding time (Day)	130	115	83	98	122	151
Storage capacity						
Racks and vial boxes	Number of racks (ea)	6/7	6/7	6/7	6/7	6/7
	Diameter of racks (mm)	82×84	82×84	142×144	142×144	142×144
	Diameter of boxes (mm)	76×76	76×76	134×134	134×134	134×134
	Number of layers per rack (ea)	5	6	4	5	6
	Number of vials (ea)	750/875	900/1050	2400/2800	3000/3500	3600/4200
25ml (single)	Number of racks (ea)	18	18	30	30	30
	Number of layers per rack (ea)	2	3	2	2	3
	Number of blood bags (ea)	36	54	60	60	90
	Number of racks (ea)	7	7	7	7	7

25ml (single)	Number of racks (ea)	18	18	30	30	30	30
	Number of layers per rack (ea)	2	3	2	2	3	4
	Number of blood bags (ea)	36	54	60	60	90	120
25ml (multi)	Number of racks (ea)	7	7	7	7	7	7
	Number of layers per rack (ea)	2	2	2	2	3	4
	Number of boxes per layer (ea)	4	4	10	10	10	10
	Number of blood bags (ea)	56	84	140	140	210	280
50ml (single)	Number of racks (ea)	—	—	—	30	30	30
	Number of layers per rack (ea)	—	—	—	3	4	4
	Number of blood bags (ea)	—	—	—	90	120	120
50ml (multi)	Number of racks (ea)	—	—	—	7	7	7
	Number of layers per rack (ea)	—	—	—	2	3	4
	Number of boxes per layer (ea)	—	—	—	8	8	8
	Number of blood bags (ea)	—	—	—	112	168	224

Containers						
Height (mm)	718	754	765	790	870	960
Outer diameter (mm)	508	461	681	681	681	681
Empty weight (kg)	15	15.2	38.3	41.2	42.3	48.9
Accessories included						
Lockable lid	√	√	√	√	√	√
Protective holding bag	—	√	—	—	—	—
Optional accessories						
Intelligent monitor	√	√	√	√	√	√
Carrying cart	√	√	√	√	√	√