## Thermo scientific

### PRODUCT INFORMATION DNA Polymerase I

Pub. No. MAN0013723 Rev. Date 03 May 2016 (B.00)

Lot: \_

#### Expiry Date: \_

Components	#EP0041 500 U	#EP0042 2500 U
Concentration	10 U/µL	10 U/µL
10X Reaction Buffer	1 mL	5 x 1 mL

#### Store at -20 °C

#### www.thermofisher.com

#### Description

DNA Polymerase I, a template-dependent DNA polymerase, catalyzes 5' $\rightarrow$ 3' synthesis of DNA. The enzyme also exhibits 3' $\rightarrow$ 5' exonuclease (proofreading) activity, 5' $\rightarrow$ 3' exonuclease activity and ribonuclease H activity.

#### Applications

- DNA labeling by nick-translation in conjunction with DNase I (1-3).
- Second-strand cDNA synthesis in conjunction with RNase H (4).

#### Source

*E.coli* cells with a cloned *polA* gene from *E.coli*.

#### **Molecular Weight**

103 kDa monomer.

#### **Definition of Activity Unit**

One unit of the enzyme catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction 30 min at 37 °C.

#### Storage Buffer

The enzyme is supplied in: 25 mM Tris-HCI (pH 7.5), 0.1 mM EDTA, 1 mM DTT and 50% (v/v) glycerol.

#### **10X Reaction Buffer**

500 mM Tris-HCl (pH 7.5 at 25 °C), 100 mM MgCl\_2, 10 mM DTT.

#### Inhibition and Inactivation

- Inhibitors: metal chelators, PP<sub>i</sub>, P<sub>i</sub> (at high concentrations) (5).
- Inactivated by heating at 75 °C for 10 min or by addition of EDTA.

#### Note

- DNA Polymerase I accepts modified nucleotides (e.g. biotin-, digoxigenin-, fluorescent-labeled nucleotides) as substrates for the DNA synthesis.
- Activity of DNA Polymerase I in Thermo Scientific buffers (in comparison to activity in assay buffer):

Buffers	Activity, %
for restriction enzymes:	
Thermo Scientific™ FastDigest™, FastDigest™	
Green, O, R, 1X Thermo Scientific <sup>™</sup> Tango <sup>™</sup> ,	
2X Tango™, BamHI, EcoRI	100
G	75-100
Ecl136II, Pacl, Sacl, Kpnl	50-75
В	25-50
for <i>Taq</i> buffer with KCl,	
Tag buffer with $(NH_4)_2SO_4$ and	100
<i>Pfu</i> buffer	
for RT buffers	100

# **CERTIFICATE OF ANALYSIS**

#### Endodeoxyribonuclease Assay

Incubation of supercoiled plasmid DNA with polymerase.

Quality authorized by:

Jurgita

Jurgita Zilinskiene

(continued on back page)

# Protocol for Radioactive DNA labelling by nick-translation

1. Mix the following components:

10X Reaction Buffer	2.5 µL	
Mixture of 3 dNTPs, 1 mM*	1.25 µL	
(without the labeled dNTP)		
$[\alpha - {}^{32}P]$ -dNTP	1.85-3.7MBq	
~110 TBq/mmol (3000 Ci/mmol)	(50-100 µCi)	
DNase I, RNase-free (#EN0521)	1	
freshly diluted to 0.002 U/µL**	1μL	
DNA Polymerase I	0.5-1.5 µL (5-15 U)	
Template DNA	0.25 µg	
Water, nuclease-free (#R0581)	to 25 µL	

- To prepare the mixture of three non-labeled dNTPs (1 mM of each), mix 1 µL aliquots of stock solutions of each dNTP (100 mM, from #R0181) with 97 µL of Water, nucleasefree (#R0581). These dNTP mixes can be stored at -20 °C for further use.
- \*\* DNase I, RNase-free (#EN0521) can be diluted with the 1X reaction buffer for DNA Polymerase I.
- 2. Immediately incubate at 15 °C for 15-60 min.
- 3. Terminate the reaction by adding 1 µL of 0.5 M EDTA, pH 8.0 (#R1021).
- Take an aliquot (1 μL) to determine efficiency of the label incorporation. A specific activity of DNA at least 10<sup>8</sup> cpm/μg DNA is expected.
- 5. If needed, the labeled DNA may be separated from the unincorporated radioactive precursors on Sephadex G-50 or Bio-Gel P-60 column.

#### Note

- The reaction volumes can be scaled up or down providing that the final concentrations of the components (DNA, dNTPs, labeled dNTP) are as indicated in the protocol.
- Radioactive DNA probes with higher specific activities can be prepared using two radioactively labeled dNTPs simultaneously. In this case, the composition of the unlabeled dNTP mix should be adjusted accordingly.

#### References

- 1. Ausubel, F.M., et al., Current Protocols in Molecular Biology, vol. 1, John Wiley & Sons, Inc., Brooklyn, New York, 3.5.3-3.5.6, 1994-2005.
- 2. Sambrook, J., Russell, D.W., Molecular Cloning: A Laboratory Manual, Third Edition, Cold Spring Horbor laboratory, Cold Spring Harbor, N. Y., 2001.
- 3. Yu, H., et al., Cyanine dye dUTP analogs for enzymatic labeling of DNA probes, Nucleic Acids Res., 22, 3226-3232, 1994.
- 4. Gubler, U., Hoffmann, B.J., A simple and very efficient method for generating cDNA libraries, Gene, 25, 263-269, 1983.
- 5. Eun, H-M., Enzymology Primer for Recombinant DNA Technology, Academic Press, INC, 1996.

# LIMITED USE LABEL LICENSE: Internal Research and Development Use Only.

The purchase of this product conveys to the buyer the limited, non-exclusive, nontransferable right (without the right to resell, repackage, or further sublicense) to use this product for internal research and development purposes. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of the product does not include or carry any right or license to use, develop, or otherwise exploit this product commercially and no rights are conveyed to the buyer to use the product or components of the product for purposes including but not limited to provision of services to a third party, generation of commercial databases or clinical diagnostics. This product is sold pursuant to authorization from Thermo Fisher Scientific and Thermo Fisher Scientific reserves all other rights. For information on purchasing a license for uses other than internal research and development purposes, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Inc., 5791 Van Allen Way, Carlsbad, California 92008.

#### PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively *for research purposes and in vitro use only.* The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to <u>www.thermofisher.com</u> for Material Safety Data Sheet of the product.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.