

**Catalog Number:** 431484

**Product Description:** Corning® Spin-X® UF Concentrator, 6 mL, 30k MWCO,  
Non-Sterile

**Component Materials:**

- Tube & Snap Lid - Virgin Polycarbonate, meets *USP, Class VI* requirements for plastic containers and closures.
- Cap - Virgin High Density Polyethylene, Heavy metal free (meets *CONEG* req.) color concentrate.
- Membrane - Polyethersulfone (PES).

**Product Dimensions:**

- Tube Height w/Cap - 4.80 inches
- Tube Diameter at Base - 0.67 inches
- Cap OD - 0.87 inches

**Maximum RCF:**

- Swing out (bucket) 4,000xg
- Fixed Angle (25°) 8,000xg

**Membrane Retention:**

The product has been tested for retention using a standard molecular weight marker solution.

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

- Integrity Test: Pressure hold test method.
- Visual Attributes: Visual examination of product.
- Packaging: Visual examination of accurate labeling and correct product configuration.

**Revision Date:**

03-14-19

Rev No: 2