

Corning® Gosselin™ Centrifuge Tubes

Designed for Optimum Performance

CORNING

Corning offers a full line of disposable, premium quality conical centrifuge tubes manufactured from high quality, ultra-clear resins. These tubes meet the quality and performance standards required by quality control laboratories.

- ▶ High speed 15 mL and 50 mL polypropylene tubes
 - 2500 to 10000 RCF rating
 - 95kPa resistant according IATA Packing Instruction 650
 - Sterility Assurance Level (SAL) of 10^{-3}
- ▶ Improved plug seal cap with leak-resistant design
- ▶ Improved packaging with resealable, easy open bulk packaged sleeves

Ordering Information

Cat. No.	Size (mL)	Sterile	Design	Cap	RCF (xg)	Qty/Pk	Qty/Cs
TC15A-03	15	No	Conical base	Plug seal	8000	25	500
TC15A-01	15	Yes	Conical base	Plug seal	8000	25	500
TC50-01	50	No	Conical base	Plug seal	10000	25	500
TC50-04	50	Yes	Conical base	Plug seal	10000	25	500
TU50-012	50	No	Self-standing	Plug seal	2500	25	500
TU50-013	50	Yes	Self-standing	Plug seal	2500	25	500



For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

CORNING

Corning Incorporated
Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

ASIA/PACIFIC
Australia/New Zealand
t 61 427286832

China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6572-9740
f 65 6735-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE
CSEurope@corning.com

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 020 655 79 28
f 020 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries
t +31 (0) 206 59 60 51
f +31 (0) 206 59 76 73

LATIN AMERICA
grupoLA@corning.com

Brasil
t 55 (11) 3089-7400
Mexico
t (52-81) 8158-8400