

Constant-Temperature Incubator

Features:

- * LCD display.
- * Microcomputer temperature controller.
- * Double door design and magnetic door sealing.
- * Castors are equipped for easy movement.
- * The working chamber adopts round angle structure, easy to clean.



Safety Functions:

- * Over temperature protection, temperature sensor malfunction alarm.
- * Parameters can be saved automatically when power off or system error.

Email: export@biobase.com

Technical Parameters:

Model	BJPX-H50II	BJPX-H60II	BJPX-H160II	BJPX-H270II
Capacity	48L	84L	160L	270L
Temp. Range	RT+5~60°C			
Temp. Fluctuation	±0.3°C			
Time Range	1~9999mins			
Shelves	2 pcs/adjustable		3 pcs/adjustable	
Consumption	437W	502W	580W	688W
Power Supply	AC220V±10%, 50/60Hz(Standard); AC110±10%, 60Hz(Optional)			
Internal Size(mm)	350*348*400	410*405*510	505*455*700	600*600*750
External Size(mm)	510*548*740	570*590*850	673*645*1033	768*795*1083
Package Size(mm)	590*850*930	730*700*1030	800*800*1210	930*940*1250
Gross Weight(kg)	50	75	99.5	150